IT

## KIT PERCUTANEO UNICO FORTY

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
### INSTRUCTIONS FOR USE

#### GENERAL DESCRIPTION OF "CENTESIS"

The procedure consists of percutaneously inserting a drainage catheter in order to draw out fluids (liquid and air) mainly but not exclusively from the pleural cavity of the patient. The purpose of this treatment is to relieve the symptomatology related to pneumothorax, including hypertensive, with pleural effusion and spillage of various nature; the procedure may also include drawing of fluids for diagnostic purposes or infusion of therapeutic solutions.

#### INTENDED USE OF THE DEVICE

The UNICO system is prevalently intended for creation of an access, in a minimally invasive way, to the pleural cavity of the patient in order to evacuate accumulated fluids and air. The device can be used for drainage by gravity or suction connected to thoracic collection systems as well as for collection of fluid samples for diagnostic purposes. The device can also be used for infusion of liquids such as saline rinsing solutions.

#### DEVICE DESCRIPTION

The device called UNICO is a percutaneous kit available in different diameters. The kit is essentially made up of the following components:

1. A main body (1) with two accesses (2) with Luer-lock connection that allow connecting drainage devices such as syringes, collection bags, etc. according to the procedures described in this instruction leaflet. A selector (3) fitted on the main body that can be turned 90° to allow drainage by gravity or injection of liquids such as rinsing solutions. The catheter (4), a 40 cm-long with depth holes and markers is securely fastened to the main body is connected to the main body, by means of a butterfly connector (5) which allows interposing an extension line (6).
2. In all the versions, the catheter has a hole for introduction of a Verres needle (7), which is used to atraumatically introduce and position the catheter. The above mentioned needle has a stylet with an atraumatic tip that protrudes from the needle end. This stylet can slide longitudinally until uncovering the cutting edge of the needle during introduction to then return to the extended position by means of a spring. The needle is equipped with a grip (8) in which a visual indicator is fitted so that the user can see the position of the atraumatic protection stylet.

#### ACCESSORIES

The following accessories are contained in the kit:

- a 2000 ml collection bag (9) fitted with an air evacuation valve (10) and an emptying valve (11);
- a syringe (12);
- a single-use scalpel (13);
- a drainage system connector (14);
- a single-use flowstop clamp (15);
- an extension line (6);
- a catheter fixing kit (16).

#### GENERAL WARNINGS AND PRECAUTIONS

- Carefully read the instructions for use before using the device.
- Use of the device is reserved for physicians and healthcare operators qualified to perform percutaneous access procedures and aware of the risks and possible implications of the above mentioned procedures.
- This device and each of its parts must be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet. The manufacturer declines all responsibility for any event deriving from improper use or different from that indicated.
- Check that the visual indicator of the protection stylet position slides smoothly and that it shows green when the stylet is in extended position. If the green indicator does not disappear during initial compression on the external surface of the thoracic wall, do not continue introducing the needle and replace the device.
- The device has been designed for use in association with Redax products and accessories.

If different products and accessories are used, the physician is responsible for checking compatibility before use.

- The device may remain in place for maximum 29 days. The duration of the drainage procedure must be determined by the physician based on clinical considerations.
- Do not use the device if its wrapping or one of its parts is damaged.
- Patients diagnosed with coagulopathy or other blood diseases must be carefully examined for the risk of haemorrhages.
- Do not in any way handle the drain or make extra holes with the aid of cutting tools.
- If securing the drain with a stitch, be extremely careful not to suture the tube or cut it with needles or sharp tools. Do not over-tighten the suture as this might cause the drain to break.
- After positioning the drain, check that the perforated section is fully inserted in the wound and that the connections are perfectly sealed.
- Manually remove the drain pulling it gently avoiding abrupt manoeuvres. Do not use metal instruments, such as forceps or devices that may cause the drain to break.
- If the drain is left in place for very long periods of time, it may be difficult to remove. Be extremely careful during the removal operations.
- Once the drain has been removed, check that it is fully intact, an essential condition to exclude that drain fragments have been left in the cavity caused by accidental rupture.
- Single-use device. Reuse may lead to alterations in performance and cross-contamination risks.

#### PREPARING FOR USE

Before using the system, carry out the preparation procedure as described below:

1. Check that the package is intact in order to ensure product functionality and sterility.
2. Open the packages contained in the kit using a sterile procedure.
3. Determine the needle insertion site, if necessary by means of X-ray or ultrasound scan.
4. Prepare a sterile field around the selected site taking the normal precautions for skin disinfection.
5. Take the main body with attached catheter out of the package;
6. Insert the needle in the catheter through the dedicated hole indicated by the arrows (Fig. 1) until the grip sits in place.

**Note:** The arrow on the grip must match the one on the catheter and the grip inclination must be directed in such a way as to follow the tube curvature (Fig. 2). The end of the needle with its protection stylet must protrude from the catheter end by a few millimetres.

7. Perform local anaesthesia of the site selected for needle introduction using the usual anaesthetics.
8. To facilitate catheter introduction, the physician may decide to make a small incision in the skin using the single-use scalpel included in the package.

#### WARNINGS

In cases of thoracentesis, the needle access point should be selected in proximity of the upper margin of the lower rib in order to avoid hitting the neurovascular bundles at the lower margin of each rib.

#### INTRODUCING THE NEEDLE AND POSITIONING THE CATHETER

After doing the preparation described above, the needle can be introduced.

**WARNING:** Although the system is fitted with numerous safety devices, this operation is not without risks and it is therefore recommended to follow the procedure described below.

1. Firmly grip the device using both hands; use the hand on the needle grip to impart the insertion thrust. The other hand should preferably be positioned in proximity of the needle

- apex to allow maintaining the direction and avoiding excessive needle insertion after passing the intercostal space.
- Rest the needle tip on the patient's skin and exercise slight pressure until the green visual indicator disappears (Fig. 1).
  - Very slowly continue with introduction until the green visual indicator reappears; this means that the end of needle has entered the pleural space and that the protection stylet has returned to its extended position performing its safety function (Fig. 2). Entry into the pleural cavity and consequent reappearance of the green indicator may be accompanied by a "click" due to the action of the spring that extends the protection stylet.
  - Orient the needle grip in such a way as to direct catheter insertion, then retract the needle and at the same time insert the catheter in the thoracic cavity (Fig. 3). The depth markers help the physician determine the positioning depth.

#### FLUID SUCTION

Having positioned the catheter, the fluid can be drawn out of the pleural cavity with the aid of a syringe and a collection bag. To connect the above mentioned accessories, open the corresponding caps on the two Luer-lock connectors on the main body. Connect the syringe to one of the two connectors and the collection bag to the other one.

**NOTE:** The physician is free to choose which connector to use for the syringe or the bag, whatever is most comfortable.

Pull the syringe plunger (Fig. 4a): the fluid will be drawn out of the pleural cavity into the syringe. Pushing on the plunger, the fluid will automatically be sent to the collection bag (Fig. 4b). Repeating the pulling and pushing operation in sequence, the fluid will be evacuated and collected in the bag.

**NOTE:** No additional operation is necessary and the selector must remain in the initial position (the knob is oriented according to the catheter axis as shown in Fig. 4).

The bag can be replaced or the syringe or both accessories removed at any time based on the drainage result. When removing the accessories, it is advisable to close the protective caps to prevent possible contamination. In any event, should the caps remain open, there is no risk of backflow to the patient.

#### FITTING THE EXTENSION LINE

The device version with butterfly connector allows fitting an extension line between the catheter and the main body. That way, the drain can easily be secured and kept in place for a long time providing greater comfort for the patient. To fit the extension line, clamp the catheter using a flowstop clamp. Then disconnect the butterfly connector and connect the line as shown in Fig. 5. After completing this operation, the clamp can be opened and you can fasten the catheter to the connector of the fluid collection accessories.

#### FLUID DRAINAGE: COLLECTION BAG

After positioning and securing the catheter, you can proceed with fluid drainage by gravity. In this case, you need to connect the drainage bag included in the package to any one of the connectors available; the drained fluids are collected in the bag by effect of gravity.

To perform fluid drainage, operate as follows:

- Connect the collection bag to one of the Luer-lock connectors and make sure that the other connector is closed with its cap.
- Appropriately position the collection system below the level of the patient.
- Turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position (Fig. 6). At this point, the catheter and the collection system are directly connected without interposition of any valve.

**WARNING:** Use only the collection bag included in the kit or Redax spare parts as they are equipped with a one-way valve that prevents fluid backflow to the patient.

- If you want to replace the collection bag, for example, when it is full, turn the selector back to the initial position and then replace the bag. Once you have replaced the bag, again turn the selector by 90°.

#### FLUID DRAINAGE: THORACIC DRAINAGE SYSTEM

After positioning and securing the catheter, you can proceed with fluid drainage by gravity or suction using a thoracic drainage

system. In this case, you need to connect the drainage system to any one of the connectors available using the dedicated connector. To perform fluid drainage, operate as follows:

- Connect the connector to one of the Luer-lock connectors.
- Prepare the drainage system following the instructions for use of the respective products and then connect the system to the free end of the connector.
- Turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position. At this point, the catheter and the collection system are directly connected without interposition of any valve.

**WARNING:** Preferably use Redax products as they are equipped with a one-way valve that prevents fluid backflow to the patient. If connecting thoracic drainage systems of different manufacturers, make sure that they are equipped with a water or other type of one-way valve in order to ensure proper evacuation from the pleural cavity.

- If you want to replace the collection system, for example, when it is full, turn the selector back to the initial position and then replace the bag or drain. Once you have completed the operation, again turn the selector by 90°.

#### CONTRAINDICATIONS

From examination of the medical literature, it results that the thoracentesis procedure may in very rare cases induce iatrogenic pneumothorax. This complication has been connected with handling errors in performing thoracentesis and using large-gauge needles and drains.

There are no specific contraindications for use of this device.

#### CLEANING AND MAINTENANCE

The device is single-use and hence does not require any cleaning or maintenance. Reuse may lead to alterations in performance and cross-contamination risks.

#### DISPOSAL

After use, the device and all its accessories may be contaminated with biological material, therefore, they must be disposed of as "biologically hazardous" waste in accordance with the regulations in force in the country where the device is disposed of.

**STERILE-Sterilized by ethylene oxide. The product is sterile if the package is intact.**

Date of issue of latest version:  
See the last page: (REV.: XX-XXXX)

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