REDAX™

IT

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INSTRUCTIONS FOR USE

GENERAL DESCRIPTION AND BACKGROUND

In cardio-thoracic surgery, it is often necessary to have a regulated suction source that, when applied to the postoperative drainage system, increases its efficiency in collecting fluids (air and fluids) from the patient's chest. In order to prevent limiting patient's ability to move around, either independently or with assistance, modern suction sources are self-contained, portable and compact devices.

DESCRIPTION AND MEANING OF SYMBOLS AND ABBREVIATIONS

cmH₂O, kPa : pressure units, available on the RedLine unit
DEVICE DESCRIPTION

The Vacuum Unit is a portable, adjustable suction system dedicated for use with compatible Redax cardio-thoracic drainage systems. The

EN

vacuum unit consists of a shockproof plastic body which houses: the user interface display (3), the keypad (2), a support surface with fastening lugs (3) to ensure drainage to the unit, the connection to the filter/valve (4), the connection for the power supply unit (5).

ACCESSORIES

The vacuum unit features the following accessories:

- a power supply unit (8) for charging the batteries to ensure continuous operation of the vacuum unit; the power supply unit is equipped with plugs (10) for different international power supply standards: please refer to the respective user manual for further details.
WARNING: use only the power supply unit supplied, the specifications of which are given in these operating instructions.
- a filter/valve (available separately ref 10539) for safe connection to the chest drainage system.

INTENDED USE OF THE DEVICE

The Dretech® REDLINE Chest Drainage device (hereinafter referred to as "Vacuum Unit") is a portable, rechargeable battery-operated unit that generates adjustable suction power. Its use is indicated following Cardio-Thoracic Surgery when a portable, self-contained vacuum source is required.

TARGET PATIENTS

The device can be used on all patients, with no restrictions on age or pathology. The attending physician must first ascertain whether there are any particular clinical situations in which the use of this device is contraindicated.

TARGET USERS

The use of this device should only be reserved for experienced medical staff who are trained to carry out chest drainage techniques and aware of potential complications associated with them. It is recommended to use the device following specific training by trained competent staff or Redax staff.

COMPATIBLE DEVICES

The device is only to be used in combination with Redax collection systems, such as Dretech Simple, Dretech Simply and Dretech Chest.

GENERAL WARNINGS AND PRECAUTIONS

- **The unit is supplied non-sterile and does not need to undergo any sterilisation cycle**
- **The vacuum unit should only be used with Redax disposable chest drainage systems for which compatibility of use is expressly indicated. For this, see the respective operating instructions.**
- **The vacuum unit must not be used if the main body and/or one of its parts is obviously damaged.**
- **Carefully follow the preparation steps before using the system.**
- **Do not expose the vacuum unit and accessories to contact with liquids or atmospheric agents in conditions that exceed the protection rating IP44. For cleaning operations, please refer to the following paragraph.**
- **This device and each of its parts may only be used in the absence of risks in the field of intended application and in the manner specified in this instruction sheet, corresponding to the product type. The manufacturer denies all liability for improper use or use other than that indicated herein.**
- **The device must only be used by qualified and specialised medical personnel (doctors and/or nurses). The manufacturer accepts no liability for use by unauthorised and unqualified personnel or for use in conditions not specified in this instruction manual.**
- **Report any serious incidents occurring with the use of the device to the manufacturer and the Competent Authority.**
- **The Vacuum Unit contains rechargeable batteries and other materials which must be disposed of and/or recycled in all countries where such processes are provided for and regulated, and in accordance with local laws.**
NOTE: to dispose of batteries properly, they must be removed from the container. Battery

replacement/removal must only be carried out by specialised technical personnel authorised by REDAX.

- **During operation, do not cover the vacuum unit and/or the charger and/or the chest drainage system with clothing, gauze, etc., to avoid overheating.**
- **The attending physician and assisting personnel should be aware of the possible implications of high suction drainage for the correct use of the Vacuum Unit with values above -30cmH₂O (-3.0 kPa), for which they are ultimately responsible.**
- **When used at high altitudes or on board aircraft, the performance of the vacuum unit may differ from the rated performance. Max. altitude : 2000m a.s.l.**
- **Do not expose to high temperatures or fire, due to the risk of battery explosion.**
- **REUSABLE DEVICE: AT THE END OF USE DO NOT THROW THE UNIT AWAY WITH THE COLLECTION SYSTEM! To prevent this, when the unit is switched off, it flashes for a certain number of times, emits an audible alarm and displays a special symbol (18) on the screen.**
- **The unit complies with all applicable electromagnetic compatibility and electrical safety requirements, in accordance with the applicable European standards (CEI EN 60601-1 series standards): in any case, care should be taken when the unit is used in the vicinity of equipment which emits electromagnetic radiation during its operation.**
- **Pay attention to the maximum temperatures that can be reached by the unit casing and power supply unit (50°C and 52°C respectively) during use.**
- **Keep the unit away from equipment used for patient treatment that may cause mutual interference.**

FIRST USE

Preliminary charging of batteries

Although the vacuum unit has a residual charge, it requires to be fully charged before it can be used for the first time.

1. Take the vacuum unit out of the packaging.
2. Plug in the socket compatible with the local standard and connect the power supply unit to a mains outlet
3. Connect the power supply unit to the vacuum unit by means of the plug located on the back of the unit (5): this will start charging the batteries in the unit.

WARNINGS for charging and batteries: the unit is equipped with rechargeable Lithium-ion batteries (Li-ion). Charging from a fully discharged battery takes approximately 2 hours.

4. Recharging is indicated by the battery symbol with an indication of the charging status in progress (*symbol d.*)

NOTE: Recharging will not take place if the batteries are damaged or have reversed polarity or are short-circuited. This will be indicated by the appropriate error signal (*symbol e.*), described in the corresponding section "Meaning of error codes".

PREPARATION FOR USE

Preparing the vacuum unit

To prepare the vacuum unit for use, proceed as follows:

1. Visually check that the unit and its accessories are intact.
2. Prepare the Redax chest drainage system for use in accordance with the respective instructions for use and connect it to the drainage catheter.
3. Assemble the Filter/Valve (8) to the RedLine and connect the unit to the drainage system as shown in Figure 1. To do this, the drainage system must be slid over the top of the vacuum unit, taking care to fit the side and rear fastening lugs of the system into the appropriate housings on the vacuum unit.

HOW TO USE THE DEVICE

Switching on the vacuum unit

Switch on the vacuum unit by pressing the I/O button (6). At power-up the unit runs a short test cycle to check the functionality of the main components. At the end of this process, the unit will start its normal operation with a suction power of 0 cmH₂O (0 kPa). If the test shows a malfunction, the corresponding error (*symbol e.*) will appear on the

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display and the unit will remain in stand-by mode. For the interpretation of error codes, see the corresponding section "Meaning of error codes" in this instruction manual.

Inactivity lock

After 10 seconds of inactivity, the unit's keypad locks to prevent accidental operation. This status is indicated by the "padlock" symbol (*symbol f.*) on the display. To unlock the keypad, press the I/O button (6) until the "padlock" symbol disappears.

Display

Under normal operating conditions, the display shows information with a level of brightness and contrast that has been designed to ensure maximum visibility in all outdoor lighting conditions. After about 1 minute of inactivity, the display enters power-saving mode, changing the character layout and the amount of information displayed. Pressing any key restores normal brightness and display mode.

Suction selection (low negativity)

The unit starts its operation with a suction level power of 0 cmH₂O (0 kPa). The suction power can be adjusted using the function keys (7). The normal low-negativity setting is in the range of -5 to -30 cmH₂O (-0.5 to -3 kPa) in steps of 5cmH₂O (0.5 kPa).

Operation by gravity

The RedLine unit is designed to provide gravity drainage functionality, without the need to be disconnected from the drainage system. If the doctor decides to suspend suction, even temporarily, it is sufficient to reduce suction to the value "0 cmH₂O", indicated on the display. In this way, the unit will automatically evacuate the air expelled by the patient and maintain the room pressure, without generating any additional suction.

At any time, the physician can restore suction by simply pressing the appropriate function keys (14).

NOTE: If the Redline is to be disconnected for traditional gravity drainage (with the collection system only), it can be switched off but **the filter/valve must be disconnected** in order to reduce pressure drops during air evacuation. Suction can be restored at any time, at the discretion of the physician, by reconnecting the filter/valve and switching the Vacuum Unit back on.

High-negativity operation

The vacuum unit normally operates within a suction range of -5 to -30cmH₂O (-0.5 to -3.0 kPa), which is universally recognised as a common range for chest applications. To increase the suction power even further, press and hold the increase button for at least three seconds until the corresponding indication on the display (symbol a.) lights up, accompanied by a short beep. At this point the lock is removed and the suction can be increased to a maximum of -60 cmH₂O (-6 kPa) in steps of 10 cmH₂O (1 kPa). Reducing the suction to -30 cmH₂O or below will automatically reset the lock.

Device shut-down

To turn the unit off, after unlocking the keypad if necessary, simply press and hold the I/O button (6) for a few seconds until the unit starts the shutdown cycle, which includes a visual and audible "Do Not Discard" warning, followed by a complete shutdown followed by the complete shutdown of the unit.

CAUTION: if the device is connected to the patient, ensure that portable suction drainage is no longer required for the patient before switching off the device. If this is not the case, the filter/valve must be removed to ensure proper continuation of gravity drainage.

Battery charging status

It is always possible to monitor the state of charge of the batteries by means of the special indication on display (*symbol c.*), which shows the charge by means of a segmented representation. The presence of all segments inside the symbol corresponds to a fully charged battery; the image of an empty battery corresponds to an almost discharged battery, so it must be recharged quickly. **NOTE:** when the battery is almost empty, a visual and audible signal is activated, indicating the need to connect the power supply. Approximately 2 hours of operation is guaranteed before the unit shuts down completely.

Periodic recharging of batteries / power supply unit.

When the batteries are empty or whenever you want to restore the full autonomy of the vacuum unit, you can recharge the batteries using the power supply unit (8). Connect the latter to the mains socket,

using the appropriate adapters according to the standards of the country where the unit is used. Then connect the appropriate plug to the vacuum unit (5). Recharging normally takes 2 hours, if batteries are completely discharged. Charging is highlighted both when the unit is running (*symbol b.*) and when the unit is switched off (*symbol d.*). It is possible to recharge the batteries during normal use of the Redline on a patient. During and after charging, the power supply unit will provide a maintenance current and supply the energy required for full functionality of the unit.

Note: remember not to place the device in a location away from the mains disconnect switch (mains socket) during charging. In this way, the isolation of the device is guaranteed.

Consumable parts: filter/valve replacement

The Filter/Valve serves to protect the Vacuum Unit from liquid ingress and contamination of the internal pneumatic circuit. Filters are supplied sterile and separate from the Vacuum Unit (for references see www.redax.it or a local Redax representative). For replacement, disconnect the filter/valve in use and dispose of it, then connect a new one. These operations must be carried out with the RedLine switched off

NOTE: the filter is disposable and for a single patient. It is strongly recommended that the Redline is replaced each time it is used on a different patient to avoid the risk of cross-contamination. In addition, the filter/valve must be replaced if necessary, e.g. if there are liquids in it or if it is used for a long time.

PROTECTIVE DEVICES

Excessive tilting

The RedLine is equipped with a tilt sensor that monitors the position of the unit during operation: if the tilt is more than 60° on any of the sides, the sensor interrupts the supply to the suction pump in order to prevent the collected liquid in the drainage from reaching the filter/valve connector and being sucked in accidentally. The unit emits an intermittent sound signal and an indication (*h. symbol*) appears on the display until the device is returned to the horizontal position. Even in this position, the pump blockage and the visual and audible signals remain active: carry out a visual check of the overall drainage system before resuming normal operation. To restore suction, press the I/O button for at least 2 seconds and the unit will return to the operating conditions prior to the inactivity lock.

WARNING: In the event of activation of the over-tilt lock, it is recommended to carefully check the condition of the chest drain to make sure that no liquid has been sucked into the filter/valve; if so, replace it before resuming suction.

Excessive negative pressure

Redline is equipped with an additional protective sensor that independently monitors the negative pressure level. In the event of a single fault or malfunction, which causes an excessive increase in suction power, the protective sensor interrupts the supply to the suction pump in order to safeguard the patient. When the protective sensor is activated, an error message (*symbol e.*) appears on the display, accompanied by an intermittent audible alarm. Normal functionality of the unit is restored when the negative pressure returns within the rated operating range. If the sensor trip occurs repeatedly, replace the unit in use and carry out a maintenance check.

Hydrophobic protection filter

The filter present in the accessory referred to as filter/valve, consists of a hydrophobic membrane capable of retaining liquids and preventing them from entering the internal pneumatic circuit. This accessory is disposable and for a single patient to avoid risks of cross-contamination. In order to keep this protective device in good working order, it must be replaced even if it is only partially invaded by liquids and after a very long period of use (more than 7 days).

Automatic positive pressure release valve

The one-way valve present in the accessory indicated as filter/valve, is designed to allow the evacuation of air in case of sudden flows that cannot be immediately evacuated by the suction pump (e.g. coughing). The valve is also a protective device against accidental switch-off of the Redline. In fact, in that case, the valve ensures that the patient's expelled air escapes, thus avoiding the risk of a pneumothorax.

Circuit leak

If there is a loss of integrity of the patient circuit or an accidental disconnection, a very high air leakage may occur, causing the suction pump to run continuously. After 2 minutes of continuous operation the RedLine activates an error signal (*symbol e.*) accompanied by an intermittent sound signal indicating that the connections and drainage

circuit should be checked immediately.

DISABLING SOUND SIGNALS: sound signals can be disabled by following these steps:

- 1) the unit must be switched off;
- 2) simultaneously press the start/stop I/O button (6) and the "UP arrow" of the suction power adjustment buttons (7), until the appropriate message (*symbol i.*) appears on the display.

From now on, all sound signals are disabled.

Note: in the case of sound deactivation, the only audible warning that remains active is that of high negativity. In all other cases, pay closer attention to the messages and indications on the display.

If you want to reactivate the sound alarms, follow these steps:

- 1) the unit must be switched off;
- 2) simultaneously press the I/O power button (6) and the "DOWN arrow" of the suction power adjustment buttons (7) until the indication on the display disappears.

MEANING OF SYMBOLS AND ERROR CODES

The error indication (*symbol e.*) indicates an abnormal operating situation that requires operator intervention

Meaning of the error codes:

ERROR 01: battery charging circuit malfunctioning

ERROR 02: inclinometer malfunctioning

ERROR 03: malfunction of protective pressure sensor

ERROR 05: error high air leak - probable disconnection/loss of patient circuit

ERROR 08: battery malfunction - cannot be recharged

ERROR 09: battery temperature too high - charging is interrupted

ERROR 12: Excessive negative pressure (protective sensor activation)

NOTE: the error codes are not consecutive as some of them are only used for internal diagnostic processes.

WHAT TO DO IF.... (TROUBLESHOOTING)

EVENT	PROBABLE CAUSE/EFFECT	ACTION
Penetration of liquids into the unit has occurred	The drainage system overturned and it was a long time before the upright position was re-established; liquid from the collection system entered the filter and then the internal circuit of the RedLine. Check for the presence of the over-tilt warning (<i>symbol h.</i>).	Replace the chest drain and filter and proceed with the extraordinary cleaning of the internal circuit, following the instructions in the dedicated section of this manual.
High air leak	Error code 05 Possible disconnection of a fitting or presence of a leak on the patient circuit.	Check all connections of the system by looking for possible disconnections or damage.
An error message (<i>symbol h</i>) appears on the display	An abnormal operation of the unit has occurred.	Check the error code in the list above
The unit won't switch on	The battery may be completely discharged.	Plug in the power supply unit and try switching the unit on again.
No sound signal can be heard	Sound signals have been disabled	Check for the presence of the appropriate symbol indicating the disabling of sound alarms (<i>symbol i.</i>). To reactivate them, follow the procedure described in the dedicated section.
Unit accidentally fell from a height > 1m	Accidental damage during normal use	Accidental dropping can cause damage that is not even visible, so it is strongly recommended to carry out an extraordinary maintenance check before putting the unit back into service.

The unit no longer drains	1. Error 12 2. Low patient leaks	1. the protective sensor has cut off the power supply to the suction pump due to excessive negative pressure. Check the drainage system for possible causes, in particular any blockage of the filter/valve pipe. If the problem recurs systematically, replace the Redline. 2. The pump does not start because the patient has very little air/liquid leaks. Check the patency of the patient tube.
Battery charging does not start or stops unexpectedly	1. Error 01 2. Error 08 3. Error 09	1. the internal circuit, which is responsible for recharging the batteries, has a malfunction and therefore does not recharge the batteries. 2. the batteries are behaving abnormally so charging has been suspended. 3. The temperature of the batteries has exceeded the permissible limit so charging is suspended. In all cases the Redline must be replaced and serviced.
Error during initial self-diagnosis test	1. Error 02 2. Error 03	These errors indicate that, during the initial test phase, an abnormal operation of the tilt sensor (error 02) or the protective pressure sensor (error 03) was detected. In this case, repeat the switch-on and if the error persists, contact the service centre.

NOTE: if the above instructions do not resolve the problem, contact technical support.

GUIDELINES FOR CLEANING THE DEVICE

For cleaning and disinfection of the external parts of the unit and its accessories, use the normal products available in the hospital, provided they are colourless and not excessively acidic or basic, as described below:

1. pour a small amount of product onto a cloth or gauze;
2. wipe the casing of the Unit in all its parts. Preferably use colourless skin disinfectants, low chlorine solutions, methyl alcohol, pH-neutral detergent solutions.

Warnings:

- a 10% dilution of sodium hypochlorite solution or similar is recommended.
- do not use disinfectant solutions and/or detergents containing dyes, in order to avoid damaging the casing or other parts of the unit and its accessories.
- do not use acidic or basic solutions other than those recommended to avoid damage to exposed parts of the unit and its accessories.
- always use cleaning and disinfecting solutions with the aid of a cloth or gauze.

NEVER POUR THE PRODUCT DIRECTLY ONTO THE UNIT AND ITS ACCESSORIES.

clean the unit after use on each patient.

EXTRAORDINARY CLEANING OF THE INTERNAL CIRCUIT

If contamination is suspected, the internal pneumatic circuit of the pump can be cleaned by flushing with white spirit (ethyl alcohol) following the steps described below:

1. Fill a syringe with 30 ml of alcohol;
2. Preferably work over a sink or washable surface

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3. Switch on the pump and set a pressure of -20cmH₂O;
4. Move the syringe closer to the filter connection site and slowly inject the alcohol;
5. The alcohol will be sucked in by the pump and will exit through the alcohol outlet at the bottom of the unit;
6. Once the suction of the liquid is complete, leave the pump on for at least 5 minutes to help dry the circuit and evaporate the alcohol completely. After 5 minutes, switch off the pump.

NOTE: perform cleaning outside the patient's environment. This operation has been defined to be carried out in exceptional cases where there has been clear contamination of the internal circuit of the RedLine and is not intended as a routine procedure, as it could cause premature deterioration of the components.

WARRANTY, MAINTENANCE AND PERIODIC CHECKS

The unit requires no routine maintenance by the user other than the cleaning operations described in the previous paragraph. All maintenance work must be carried out by authorised Redax technicians.

In the event of obvious damage to the display or the casing, or if the battery charge duration is insufficient, replace the device and contact the authorised Redax network for maintenance.

STORAGE AND CONSERVATION

It is recommended that the unit be stored in its own case to prevent damage from accidental drops. From the point of view of chemical and physical stability, the product and the material from which it is made do not change over time. Exposure to sunlight or artificial light sources does not change the structure of the product. However, storage at room temperature is recommended and exposure to high temperatures and ultraviolet radiation should be avoided.

DISPOSAL

The Unit contains electronic components and rechargeable lithium-ion batteries, as well as the materials making up the casing and other supporting parts. For this reason, the device must be disposed of in accordance with the applicable regulations for the disposal of electronic waste.

For the European Community

Drehtech® Redline is marketed in compliance with the WEEE (Waste Electrical and Electronic Equipment) regulation.

The class II device has no applied parts but has been subjected to all tests related to the TYPE B APPLIED PART in accordance with CEI EN 60601-1.

CONTRAINDICATIONS:

There are no known contraindications to using the RedLine device for its intended purpose.

Symbols used

The unit has been labelled in accordance with current European standards, in particular according to CEI EN 60601-1 and related standards. In addition, symbols have been used which are not specified by international directives or standards, but are universally recognised or described at the beginning of this manual.

SYMBOL	DESCRIPTION
	Class II device (Ref. IEC 60417-5172)
	Type B applied part (UNIT CASING) (Ref. IEC 60417-5840)
	Refer to the instruction booklet (Ref. ISO 7010-M002)
	Separate disposal of electrical and electronic components (Ref. 2002/96/EC and 2006/66/EC)
	Direct current (Ref. IEC 60417-5031)

	Enclosure protection rating (Ref. CEI EN 60529): - protected against foreign bodies with a diameter of 1 mm or more - protected against splashing water
	Date of production (Ref. ISO 15223-1)
	Manufacturer (Ref. ISO 15223-1)
	Serial number (Ref. ISO 15223-1)
	Catalogue number (device code) (Ref. ISO 15223-1)
	Consult instructions for use (Ref. ISO 15223-1)
	Medical device (Ref. ISO 15223-1)

In addition, these symbols, which are not specified by international directives or standards but are universally recognised, have been used:

SYMBOL	DESCRIPTION
	Universal recycling symbol
	Unofficial, but universally recognised symbol referring to the European Directive 2011/65/EU

Below are the transport and storage symbols on the unit's outer box.

SYMBOL	DESCRIPTION
	Latex free
	Do not use if packaging is damaged
	Do not expose to direct sunlight
	Fragile, handle with care
	Keep dry
	Temperature limits 10 to 35 °C
	Upper side
	Medical device; US federal law restricts this device to sale by or on the order of a physician.

TECHNICAL FEATURES

VACUUM UNIT

Power supply: battery pack 11.1 V 3s1p Li-Ion 3500 mAh
Operating temperature: 15°C – 35°C
Storage temperature: -20°C - +30°C
Humidity (operation/storage): 30% - 70%
Atmospheric pressure (operation/storage): 70.0 – 106.0 kPa
Measurement tolerance: suction: ± 10%
Weight: 800 g
Protection rating Ip44

POWER SUPPLY UNIT 8A

Model: Powerbox EXM30 5009
Power supply (Volt): 15VDC
Operating temperature: 0°C – 40°C
Storage temperature: -40°C - +75°C
Max. power consumption: 30W
Maximum output current: 2.45A

POWER SUPPLY UNIT 8B

Model: TR30RDM150
Power supply (Volt): 15VDC
Operating temperature: 0°C – 40°C
Storage temperature: -40°C - +75°C
Max. power consumption: 30W
Maximum output current: 2.4A

KEY

- 01 Chest drainage fastening point
- 02 Keyboard
- 03 Display
- 04 Connection for filter/valve
- 05 Connection for power supply unit
- 06 I/O key
- 07 Suction power adjustment keys
- 08 A - B: Power supply unit (Accessory)
- 09 Power cord
- 10 Electrical plug (available in different versions)

ANNEX: Non-harmonised symbols used in the RedLine

SYMBOL	DESCRIPTION
	high suction power operation symbol
	Recharging in progress / power supply connected symbol with unit switched on
	Battery status bar symbol
	Charging in progress symbol / power supply connected with unit switched off
	error symbol
	Padlock symbol (keypad lock)
	DO NOT DISPOSE symbol
	over-tilt symbol
	sound alarms disabled symbol

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