Figura 3

REDAX™

IT

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INSTRUCTIONS FOR USE

DECLARATION OF LIABILITY

This document is the exclusive property of Redax S.p.A. All rights reserved. Any disclosure, reproduction or transfer of content to third parties is prohibited without the prior authorization of the Company. Redax S.p.A. denies all liability with respect to damage to people or property due to improper use of this product and failure to observe the indications, warnings, instructions and precautions in this user manual.

In the event of a serious accident related to the medical device, the user must notify the manufacturer and the competent authority of the state in which the event occurred.

This user manual is supplied in printed form and accompanies the device.

MANUFACTURER

Redax S.p.A.
Via Galileo Galilei, 18, 46025 Poggio Rusco (MN) Italia
P.IVA: 02556750368

info@redax.it
www.redax.it

DESCRIPTION AND MEANING OF SYMBOLS AND ABBREVIATIONS

- **cmH2O, kPa:** pressure unit (possible user choice)
- **ml:** unit of measurement of the volume of drained liquid
- **ml/min:** unit of measurement of the flow of air from the patient's chest
- **min (minutes), h (hours):** units of time used in the measurement of clinical parameters

1. DEVICE DESCRIPTION

Drentech® SilverLine (part number 10198) is a digital, transportable, rechargeable battery-operated chest suction system capable of autonomously generating adjustable suction from 0 to -60cmH2O, independent of centralised vacuum and/or power sources used following thoracic surgery and cardiac surgery.

Drentech® SilverLine can be used with compatible Redax drainage systems such as Drentech® Simple Silver, Simply and Simply Mini.

The unit allows clinical drainage parameters such as the patient's air and liquid leaks and intrapleural pressure to be displayed on a digital display Two additional functionalities have been integrated into the device: wireless connection via Bluetooth technology to a mobile device through a dedicated application allows the transmission, display and storage of recorded data.

2. INDICATIONS FOR USE

The Drentech® SilverLine thoracic suction device is a digital drainage system capable of draining fluids (air or liquid) through suction in post-operative thoracic surgery applications.

PATIENT GROUP

The device can be used on all patients, without age limitation.

USER GROUP

The device must be used by qualified medical personnel, such as surgeons and nurses.

CONTRAINDICATIONS

There are no known contraindications to using the SilverLine device for its intended purpose.

WARNINGS

- The unit is supplied non-sterile and does not need to undergo any sterilisation cycle.
- The vacuum unit must only be used with Redax-compatible disposable chest drainage systems. For a list of compatible devices, see the relevant chapter.
- The vacuum unit must not be used if the main body and/or one of its parts is obviously damaged.
- No modifications to the device are permitted.
- Carefully follow the preparation steps before using the system.
- Do not expose the vacuum unit and accessories to contact with liquids or atmospheric agents in conditions that exceed IP22D. For cleaning operations, please refer to the following paragraph.
- This device and each of its parts may only be used in the field of application and in the manner specified in this instruction sheet, corresponding to the product type. The manufacturer denies all liability for improper use or use other than that indicated herein.
- The device must only be used by qualified and specialised medical personnel (doctors and/or nurses). The manufacturer accepts no liability for use by unauthorised and unqualified personnel or for use in conditions not specified in this instruction manual.
- However, the user must always consider the drainage status indications provided by the collection system to assess the clinical condition of the patient during the drainage procedure.
- Report any serious incidents occurring with the use of the device to the manufacturer and the Competent Authority.
- The Vacuum Unit contains rechargeable batteries and other materials which must be disposed of and/or recycled in all countries where such processes are provided for and regulated, and in accordance with local laws.
- To recharge the vacuum unit battery, use the charger supplied by the manufacturer.
- To dispose of batteries properly, they must be taken out of the container. Battery replacement/removal must only be carried out by specialised technical personnel authorised by REDAX.
- During operation, do not cover the vacuum unit and/or the charger and/or the chest drainage system with clothing, gauze, etc., to avoid overheating.
- The responsible physician and assistant personnel must be aware of the possible implications of high suction drainage for the correct use of the vacuum unit with values above -30cmH2O (-3.0 kPa), for which they are ultimately responsible.
- When used at high altitudes or on board aircraft, the performance of the vacuum unit may differ from the standard. Max. altitude: 2000m a.s.l.
- Do not expose to high temperatures or fire, due to the risk of battery explosion.
- REUSABLE DEVICE: AT THE END OF USE DO NOT THROW THE UNIT AWAY WITH THE COLLECTION SYSTEM! To

prevent this, the unit emits an audible alarm when it is switched off and displays an appropriate symbol(symbol g) on the screen.

- The unit complies with all applicable electromagnetic compatibility and electrical safety requirements, in accordance with the applicable European standards (CEI EN 60601-1 series standards); in any case, care should be taken when the unit is used in the vicinity of equipment which emits electromagnetic radiation during its operation.
- Pay attention to the maximum temperatures that can be reached by the unit casing and power supply unit (50 °C and 52 °C respectively) during use.
- Keep the unit away from equipment used for patient treatment that may cause mutual interference.
- Exclusively in the event of a screen lock (03), which makes it impossible to change the negative pressure or to display the various parameter adjustment pages via the keypad (keys 08 - 09), resulting in inactivity of the suction pump, press and hold down the on/off button and unlock screen (07) for 5 seconds. This operation restores the correct functioning of the device, guaranteeing the continuity and quality of treatment.
- During the initial pairing procedure between unit and App, when prompted, enter the code **200009**

3. CLINICAL BENEFIT

The clinical benefit associated with Drentech® SilverLine is indirect: the device has an effect on the drainage procedure. The main clinical benefits of Drentech® SilverLine are:

- a) The ability to transport the unit improves patient management by facilitating the operations of hospital staff during routine activities in the hospital environment;
- b) The digital structure of the drainage system helps to reduce the duration of the drainage procedure, with a reduction in chest hose length and hospital stay compared to traditional drainage systems.

4. INSTRUCTIONS FOR USE

PRE-CHARGING OF BATTERIES - FIRST USE

The device requires to be fully charged before it can be used for the first time, despite having a residual charge. Below are the steps to be followed when using it for the first time.

1. Take the vacuum unit out of the packaging.
2. Connect the power supply to the power socket, using power adapters if necessary. The AC adapter plug is considered the means of disconnecting the device from the power supply.
3. Connect the power supply unit to the vacuum unit by means of the plug located on the back of the unit (06): this will start charging the batteries in the unit.
4. Recharging is indicated by the battery symbol with an indication of the charging status in progress (symbol d)
5. When charging is complete, replace the white protective cap.

The estimated time for a full charge of fully discharged batteries is about 2 hours.

If batteries are damaged, have reversed polarity or are short-circuited, recharging will not be possible. The device will produce the specific error signal (symbol e). Details of the error can be found in the "Meaning of Error Codes" section of this document.

PREPARATION FOR USE - CONNECTING AND SETTING UP THE SYSTEM

To prepare the vacuum unit for use, proceed as follows:

1. Visually check that the unit and its accessories are intact.
2. Prepare the Redax chest drainage system for use in accordance with the respective instructions for use and connect it to the drainage catheter.
3. Assemble the filter-(04) to the unit and connect the unit to the drainage system as shown in Figure 2 (14). The drainage system must be placed on the Drentech® SilverLine by sliding it over the top of the vacuum unit. Correct positioning is achieved by fitting the drainage system into the side and rear mounting tabs on the vacuum unit.
4. After opening the protective cap of the male luer lock connection of the vacuum unit (05), connect the luer lock connection to the end of the two-lumen hose as shown in Figure 2 (15).
5. Before activating device operation, the type of disposable system must be selected (see Chapter 6 'Compatible Devices'). To do so, proceed as follows:
 - a. Simultaneously press the I/O button (07) and the screen change button (09) until the compatible disposables to be used appear (Simple Silver, Simply and Simply mini)
 - b. Use the selection buttons (number 08) to select the system used

- c. Press and hold the I/O button (07) to confirm your choice, until the unit automatically switches off and on.

HOW TO USE THE DEVICE

STARTING THE VACUUM UNIT

To switch on Drentech® SilverLine press the I/O button (07). During start-up, the unit runs a short test cycle to check the functionality of the main components. Once the test is completed, the unit will start its normal operation with suction of 0 cmH2O (0 kPa). If the test detects system faults, the unit will remain in stand-by mode and the error symbol (symbol e) will appear on the display. For details on error codes, see the corresponding section "Meaning of Error Codes" in this operating manual.

NEW PATIENT' SELECTION (FIGURE 3-A)

When starting up the device, the display shows the 'patient' symbol (symbol l), giving you the choice if you are starting a procedure on a new patient.

The 'new patient' information in the unit can be selected within 10 seconds by selecting 'YES' using the selection buttons (08). It is not necessary to confirm the choice, the selection cursor must be left on the 'YES' indication.

The system will automatically select 'NO' if no key is pressed within 10 seconds after start-up, and will continue recording the parameters of the last recorded patient.

SUCTION SELECTION (LOW NEGATIVITY) (FIGURE 3-B)

The unit starts its operation with a suction level of 0 cmH2O (0 kPa). The suction can be adjusted using the selection buttons (08). The normal low-negativity setting is in the range of -5 to -30 cmH2O (-0.5 to -3 kPa) in steps of 5cmH2O (0.5 kPa). The unit is set to display data in cmH2O. It is possible to change the units of measurement, instructions can be found in the reference section of this operating manual.

OPERATION BY GRAVITY

Drentech® SilverLine is designed to provide drop-in drainage functionality without the need to be disconnected from the drainage system. If the doctor decides to suspend suction, even temporarily, simply reduce the suction to the value '0 cmH2O', indicated on the display.

With this setting, Drentech® SilverLine will independently evacuate the air expelled by the patient and maintain ambient pressure, without creating any additional suction.

At any time, the physician can restore suction by simply pressing the appropriate function keys.

It is possible to switch to the traditional gravity drainage system (collection system only) by switching off Drentech® SilverLine disconnecting the filter and removing the unit, in order to reduce pressure leaks in air evacuation.

At the doctor's discretion, suction can be restored at any time by reinserting Drentech® SilverLine, reconnecting the filter and switching the unit on again.

HIGH-NEGATIVITY OPERATION

The vacuum unit normally operates within a vacuum range of -5 to -30cmH2O (-0.5 to -3.0 kPa): To increase the suction intensity, press and hold the increment button for three seconds until the corresponding visual indication on the display (symbol a) and audible 'beep' sound come on. At this point, the suction can be increased to a maximum of -60 cmH2O (-6 kPa) in steps of 10 cmH2O (1 kPa).

Reducing suction to -30 cmH2O or lower will automatically reset the inactivity lock.

UNIT SELECTION

The pressure parameter can be displayed in two units: cmH2O or kPa. To change the unit of measurement, follow the steps below:

1. The unit must be switched off;
2. Press the I/O button (07) and both selection buttons (08) simultaneously to access the settings menu,
3. The settings menu presents several parameters with progressive numbering P0, P1, P2,... Scroll through the parameters by simultaneously pressing the I/O button (07) and one of the selection buttons (08), selecting parameter P1.
4. Using the selection keys (08), make your choice of unit (cmH2O, kPa).
5. To confirm your choice, press the I/O button (07) for at least 3 seconds. The unit will then switch off and start up again in operating mode, adopting the chosen pressure unit.

INACTIVITY LOCK

After 10 seconds of inactivity, Drentech® SilverLine blocks the ability to change the negative pressure via the selection buttons (08) on the main screen to prevent accidental operation. The lock status is indicated by the 'lock' symbol (symbol f) on the main screen. To release the inactivity lock, press the I/O button (07) until the 'lock' symbol disappears.

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DISPLAY

When operating Drentech® SilverLine, the information shown on the display appears with a level of brightness and contrast designed to allow maximum visibility in all outdoor lighting conditions. After about 1 minute of inactivity, the main screen enters power-saving mode by changing the character layout. By pressing any key, the brightness and graphic display are restored.

CHANGE DISPLAY SCREEN

Pressing the button at the top left of the display (09) changes the screen display to the drainage parameter display for monitoring the procedure in the patient. By pressing the (09) button once, it is possible to switch from the main screen to the sequential display of the following information:

- (1) patient's air leakage,
- (2) collection of drained fluid,
- (3) collection of fluid drained in the last 24 hours,
- (4) intrapleural pressure.

A further press of button (09) returns you to the start screen. The optional screens containing graphs of drained fluid collection and air leakage can be displayed.

AIR LEAKAGE (FIGURE 3 - C)

In this screen there is a multiple display of data.

- (1) The top shows the counting of treatment hours since the unit was switched on, reset and connected to the patient.
 - (2) The first line displays the updated air leak. The value refers to the volume of air measured by the unit in the last elapsed minute and is expressed in ml/min. The figure is updated every minute.
 - (3) In the line below, the average amount of air detected in the last hour is displayed, identified by 'AVG - 1H'. The figure is updated every hour.
- It is possible to see the leaks of the previous hours by pressing the corresponding selection buttons on the right side of the display (08): button '-' scrolls the hours backwards; button '+' returns to the current time.

DRAINED LIQUID COLLECTION (FIGURE 3 - D)

In this screen, a multiple display of data appears on the display:

- (1) The top shows the counting of treatment hours since the unit was switched on, reset and connected to the patient.
 - (2) The first line displays the level of fluid collected in the chest drain currently connected to the SilverLine. The volume of liquid is expressed in ml.
 - (3) The line below displays the total volume of liquid collected since the start of the treatment, identified by 'TOT'. The total collection of drained fluid refers to the sum of any collection systems replaced during the course of treatment, which is why it may differ significantly from the figure shown in the top line.
- WARNING:** In order to correctly detect the level of liquid collected, the unit must be placed on a flat surface. At inclinations greater than 7°, the liquid value shown on the display will be constant and will not be updated until the unit is correctly positioned on the base.

LIQUID COLLECTION DRAINED IN THE LAST 24 H (FIGURE 3 F)

In this screen there is a multiple display of data:

- (1) Above is a count of the hours of treatment since the unit was switched on, reset and connected to the patient.
- (2) In the centre is the value of the fluid drained by the patient in the last 24 hours of monitoring, expressed in ml. The figure is not available until 24 hours after the start of treatment and is updated every hour from the 25th hour.
- (3) In the line below, the words 'LAST 24H' are shown, identifying the parameter displayed.

INTRAPLEURAL PRESSURE (FIGURE 3 - E)

In this screen you can view the patient's intrapleural pressure. The parameter is detected by the pressure sensor in the unit.

- (1) On the left is an indicator for the qualitative display of intrapleural pressure: the detection range is -30 / +15 cmH2O (-3 + 1,5 kPa). Through the indicator it is possible to observe the pressure fluctuations during the patient's respiratory acts.
- (2) On the right is the graph representing the intrapleural pressure trend over the last few hours. The two lines in the graph represent the trends in maximum and minimum pressures measured during the patient's respiratory acts. You can view the graph of the previous (-) or next (+) hours by moving with the buttons (08).

NOTE: for information on the resolution data of the measurements taken by the unit, see the "Technical Specifications" section at the end of this instruction manual.

SETTINGS MENU FUNCTIONALITY

DISPLAYING OPTIONAL SCREENS

The display of Drentech® SilverLine can show the two graphs of liquid leakage and air leakage by activating the optional display mode.

The optional display mode can be activated by performing the following steps:

1. The unit must be switched off
2. Press the I/O button and both selection buttons (8) simultaneously until the settings menu appears.
3. The settings menu presents several parameters with progressive numbering P0, P1, P2... Scroll through the parameters by simultaneously pressing the I/O button (07) and one of the selection buttons (08), selecting parameter P9.
4. Using the selection buttons (8), make the choice of displaying the optional graphics. Select 'Yes' to display the graphs; confirm your choice by leaving the cursor on the 'YES' indication.
5. Press the I/O button (07) for at least 3 seconds. The unit will then switch off and start in operating mode starting with the selection of the new patient.
6. Pressing the display screen change button (09) will switch to the display of three different graphs in succession:
 - a. Qualitative trend of the sum of the patient's liquid leaks, with the possibility of also seeing the trend of the previous hours by pressing the selection buttons (8): "-" button scrolls backwards through the hours, "+" button returns to the current hour.
 - b. Qualitative trend of the patient's liquid leaks hour by hour, with the possibility of also seeing the trend of the previous hours by pressing the selection buttons (8): "-" button scrolls backwards through the hours, "+" button returns to the current hour.
 - c. Qualitative trend of the patient's air leakage with the possibility of also seeing the trend of the previous hours by pressing the selection buttons (8): "-" button scrolls backwards through the hours, "+" button returns to the current time.

DEVICE SHUT-DOWN

To switch off the device, press and hold the I/O button (6) for a few seconds. In case the keyboard is locked, proceed with unlocking first. The unit will start the shutdown cycle by a visual 'Do not throw' symbol (symbol g) and an audible signal, followed by the complete shutdown of the unit.

NOTE: In the case of patient connection, ensure that portable suction drainage is no longer required before switching off. If drainage treatment is not to be stopped but the unit is to be switched off, the filter must be removed to ensure proper continuation of gravity drainage.

BATTERY MANAGEMENT

BATTERY CHARGING STATUS

The charge status of the batteries is indicated by means of a segment representation and can be monitored via the charge indicator (symbol c). Full charge of the device battery is indicated by the presence of all segments in the charge indicator. A low battery indication is indicated by the absence of segments in the charge indicator; in this case it will be necessary to charge the device quickly to prevent the unit from shutting down.

NOTE: when the battery is almost empty, a visual and audible signal is activated, indicating the need to connect the power supply. From the time of the signal, operation is guaranteed for about 2 hours before the unit shuts down completely.

PERIODIC RECHARGING OF BATTERIES / POWER SUPPLY UNIT.

To recharge the batteries, connect the power supply unit (11) to the power socket, using if necessary the power socket adapters provided in the country in which the device is used. Connect the plug to the vacuum unit (5). Recharging normally takes 2 hours with a fully discharged battery. The charging status is highlighted with two different symbols if the device is in operation (symbol b) or not (symbol d). Battery recharging is also possible during operation of the Drentech SilverLine during on-patient drainage treatment. During and after charging, the power supply unit will provide a maintenance current and supply the energy required for full functionality of the unit.

CONSUMABLE PARTS

REPLACEMENT OF THE DISPOSABLE COLLECTION SYSTEM

When you want to replace the disposable collection system, proceed as follows:

1. Set suction value to 0 cmH2O
2. Clamping the patient's drainage hose
3. Disconnect the filter from the collection system
4. Disconnect the patient hose from the collection system
5. Remove the exhausted system, use a new collection system as a

- replacement, choosing from compatible ones
- Restore connections to the filter and patient hose
 - Remove the clamp from the patient's drainage hose
- NOTE:** The current liquid leak value resets to zero each time the collection system is replaced, while the "TOT" value is updated by adding up the amount of liquid previously recorded.

FILTER REPLACEMENT

The filter protects Drentech SilverLine from liquid ingress and contamination of the internal pneumatic circuit. For replacement, switch off Drentech SilverLine, disconnect the filter in use and dispose of it. Then connect the new filter.

Filters are supplied sterile and separate from the Vacuum Unit (for references see www.redax.it or a local Redax representative).

NOTE: The filter is disposable and its use is single-patient. It is strongly recommended that the filter be changed each time Drentech SilverLine is used on a different patient. This avoids the risk of cross-contamination. The filter must also be replaced if necessary (e.g. presence of liquids inside, prolonged use).

PROTECTION SYSTEMS EXCESSIVE TILTING

Drentech SilverLine is equipped with a tilt sensor that monitors the position of the unit during operation. If the inclination is more than 60° on any one side, the sensor interrupts the suction pump supply to prevent accidental suction of the liquid collected in the drainage.

In the event of excessive tilting, the unit emits an intermittent beep and visual signal on the display (h symbol) until the horizontal position of the unit is restored. Once the horizontal position has been restored, the pump lock and the visual and audible signals remain active: carry out a visual check of the drainage system as a whole before resuming normal operation. Press the I/O button when the indicator on the display is present to remove the inactivity lock and restore suction.

ATTENTION: In the event of a blockage due to excessive tilt, it is recommended to carefully check the condition of the chest drainage. In particular, check that no liquid has been sucked into the filter; if so, replace it before restarting suction.

EXCESSIVE NEGATIVE PRESSURE

Drentech SilverLine is equipped with an additional protective sensor that independently monitors the negative pressure level. In the event of a single fault or malfunction, which causes an excessive increase in suction, the protection sensor interrupts the suction pump supply in order to safeguard the patient. When the protective sensor is activated, an error message (symbol m) appears on the display, accompanied by an intermittent audible alarm and the word PRESS. Normal functionality of the unit is restored when the negative pressure returns within the rated operating range.

If the sensor trip occurs repeatedly, replace the unit in use and carry out a maintenance check.

HYDROPHOBIC PROTECTION FILTER AND AUTOMATIC POSITIVE PRESSURE RELEASE VALVE

The filter is a compatible accessory (see Chapter 5 "Compatible Accessories") It consists of two main parts:

- Hydrophobic membrane: retains liquids and prevents them from penetrating the internal pneumatic circuit.
- One-way valve: allows the evacuation of excess air. In the event of sudden flows that cannot be immediately evacuated by the suction pump (e.g. coughing fit), the valve serves as a protective device. It prevents problems with accidental switch-off of Drentech SilverLine.

NOTES: The valve should by no means be understood as a gravity drainage valve, it is a protective device for the device and its internal parts.

CIRCUIT LEAK

In the event of loss of integrity of the patient circuit or accidental disconnection, continuous operation of the suction pump may occur due to a very high air leakage, which is not related to a consistent air leakage of the patient. Drentech SilverLine activates an error signal visually via "LEAK" (m symbol) and audibly via intermittent sound to indicate continuous operation for 2 minutes. In this case, the connections and the drainage circuit must be checked immediately.

DISABLING SOUND ALARMS

The device allows sound alarms to be disabled. This functionality can be activated by following these steps:

- The unit must be switched off;
- Simultaneously press the ignition I/O button (07) and the "UP arrow" of the suction adjustment buttons (08) until the indication (symbol i) appears on the display.

Warning: The high negative pressure use signal (< of - 30

cmH2O), the low battery signal and the acoustic tipping signal remain active. In all other cases, pay closer attention to the messages and indications on the display.

To reactivate sound alarms, follow the steps below:

- The unit must not be in operation;
 - Simultaneously press the ignition I/O button (07) and the "DOWN arrow" of the suction adjustment buttons (08), until the indication (symbol i) on the display disappears.
- From this moment on, all sound signals are reactivated.

MEANING OF SYMBOLS AND ERROR CODES

The error indication (symbol e) indicates an abnormal operating situation that requires operator intervention Below are the error codes with their meanings.

Meaning of the error codes:

ERROR 01	battery charging circuit malfunction
ERROR 02	inclinometer malfunction
ERROR 03	malfunction of protective pressure sensor
ERROR 08	battery malfunction - unable to recharge
ERROR 09	excessive battery temperature - charging is interrupted

NOTE : the error codes are not consecutive as some of them are only used for internal diagnostic processes.

WHAT TO DO IF.... (TROUBLESHOOTING) EVENT	PROBABLE CAUSE/ EFFECT	ACTION
Penetration of liquids into the unit has occurred	Possible overturning of the system and remaining in a non-vertical position. Liquid ingress into the filter and internal circuit of Drentech SilverLine. Check for the presence of the over-tilt warning (symbol h).	Replacement of chest drain and filter.
High air leak LEAK lettering and the m symbol on the display	Possible disconnection of a connection or presence of a leak on the patient circuit.	Check all system connections for possible disconnections or damage.
Error message (symbol e) on display	An abnormal operation of the unit has occurred.	Check the error code in the list above.
The unit won't switch on.	The battery may be completely discharged.	Connect the power supply and try to switch on the unit
No sound signal can be heard	Sound signals have been disabled	Verify the presence on the display of the appropriate symbol indicating that audible alarms are disabled (symbol i). To reactivate them, follow the procedure described in the section "Disabling sound alarms".
The unit accidentally fell from a height > 1m	Accidental damage during normal use	Accidental dropping can cause even non-visible damage. An extraordinary maintenance check is strongly recommended before putting the unit back into service.
The unit no longer drains	1. The display shows the word PRESS and the symbol m	1. The protection sensor cut off the power supply to the suction pump due to excessive negative pressure. Check the drainage system for blockages in the filter hose. If the problem recurs systematically, replace Drentech SilverLine.
	2. Low patient leaks	2. The pump does not start because the patient has very little air/liquid leaks. Check the patency of the patient hose.

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WHAT TO DO IF.... (TROUBLESHOOTING) EVENT	PROBABLE CAUSE/EFFECT	ACTION
The liquid leakage display screen appears flashing and with an audible signal	The drainage collection system is almost full	Check the filling of the drainage system and replace it with a new one if necessary. The alarm will disappear automatically with the replacement of the new system.
Constant intrapleural pressure for at least 12 minutes: the intrapleural pressure display screen flashes with a sound signal	Possible blockage/obstruction in the patient hose and/or drainage catheter	Check for clots and/or other causes of obstruction in the patient hose and/or chest catheter. Check the connection between the end of the two-lumen hose and the male luer lock connection of the vacuum unit.
Battery charging does not start or stops unexpectedly: the display shows (symbol e) with an error code.	1. Error 01: Internal battery charging circuit malfunctions and does not recharge 2. Error 08: Batteries have abnormal behaviour, so charging has been suspended 3. Error 09: Battery temperature has exceeded the permitted limit: charging is suspended	Drentech SilverLine must be replaced and serviced in all cases.
Error during initial self-diagnosis test	1. Error 02: Tilt sensor error detected during initial phase tests. 2. Error 03: Protective pressure sensor error detected during initial phase tests	Repeat ignition. If the error occurs again, it will be necessary to contact Redax support.

NOTE: If the above situations do not allow the problem to be solved, always contact technical support.

5. COMPATIBLE ACCESSORIES

Hydrophobic protection filter. The filter is single-use and single-patient to avoid risks of cross-contamination. In order to maintain its efficiency, the filter must also be replaced if it is partially invaded by liquids and after very prolonged use (more than 7 days).

6. COMPATIBLE DEVICES

Drentech® SilverLine can be used with compatible Redax drainage systems:

- Drentech® Simple Silver
- Drentech® Simply
- Drentech® Simply Mini

7. GUIDELINES FOR CLEANING THE DEVICE

The cleaning procedure of the unit must be performed at the end of use on each patient, with the device disconnected from the mains. Cleaning must be carried out by gently wiping the surface of the device using a disinfectant wipe with 70% isopropyl alcohol. After cleaning, the device and its components must be inspected for signs of deterioration that may limit the lifespan and/or performance of the device, such as visible corrosion, mechanical wear, abrasion, damage or deformation. Dispose of the device if there are obvious signs of deterioration.

8. WARRANTY, MAINTENANCE AND PERIODIC CHECKS

The unit requires no routine maintenance by the user other than the cleaning operations described in the previous paragraph. All maintenance work must be carried out by authorised Redax technicians.

In the event of obvious damage to the display or the casing, or if the battery charge duration is insufficient, replace the device and

contact the authorised Redax network for maintenance.

9. STORAGE

It is recommended that the unit be stored in its own case to prevent damage from accidental falling. However, storage at room temperature is recommended and exposure to high temperatures and ultraviolet radiation should be avoided.

10. DISPOSAL

The unit contains electronic components and rechargeable lithium-ion batteries.

For this reason, the device must be disposed of in accordance with the applicable regulations for the disposal of electronic waste.

11. LABELLING AND SYMBOLS USED

The unit has been labelled in accordance with current European standards, in particular according to CEI EN 60601-1 and related standards. In addition, symbols have been used which are not specified by international directives or standards, but are universally recognised or described at the beginning of this instruction manual.

Below are the transport and storage symbols on the unit's outer box.

SYMBOL	DESCRIPTION
	Class II device (Ref. IEC 60417-5172)
	Type B applied part (UNIT CASING) (Ref. IEC 60417-5840)
	Refer to the instruction booklet (Ref. ISO 7010-M002)
	Separate disposal of electrical and electronic components (Ref. 2002/96/EC and 2006/66/EC)
	Direct current (Ref. IEC 60417-5031)
	Enclosure protection rating (Ref. CEI EN 60529): - protected against foreign bodies with a diameter of 1 mm or more - protected against splashing water
	Date of manufacture (Ref. ISO 15223-1)
	Manufacturer (Ref. ISO 15223-1)
	Serial number (Ref. ISO 15223-1)
	Catalogue number (device code) (Ref. ISO 15223-1)
	See instructions for use (Ref. ISO 15223-1)
	Medical device (Ref. ISO 15223-1)

In addition, these symbols, which are not specified by international directives or standards but are universally recognised, have been used:

SYMBOL	DESCRIPTION
	Universal recycling symbol
	Unofficial, but universally recognised symbol referring to the European Directive 2011/65/EU

Below are the transport and storage symbols on the unit's outer box.

SYMBOL	DESCRIPTION
	Latex free
	Do not use if packaging is damaged
	Do not expose to direct sunlight
	Fragile, handle with care
	Keep dry
	Temperature limits 10 to 35 °C
	Top side
	Medical device; US federal law restricts this device to sale by or on the order of a physician.

12. TECHNICAL SPECIFICATIONS

VACUUM UNIT

Power supply: battery pack 11.1 V 3s 1p Li-Ion 3500 mAh
Operating temperature: 10°C – 35°C
Storage temperature: -20°C - +30°C
Humidity (operation/storage): 30% - 70%
Atmospheric pressure (operation/storage): 70.0 – 106.0 kPa
Negative pressure measurement accuracy: ± 2.5 cmH₂O
Intra-pleural pressure measurement accuracy: ± 0.5 cmH₂O
Air leak measurement accuracy: ± 50 ml/min
Liquid level measurement accuracy:
Drentech Simple Silver: ± 15 ml
Drentech Simply: ± 10 ml
Drentech Simply mini: ± 5 ml

Weight: 855 g
Protection rating: IP22D

POWER SUPPLY UNIT 11A

Model: Powerbox EXM30 5009
Power supply (Volt): 15VDC
Operating temperature: 0°C – 40°C
Storage temperature: -40°C - +75°C
Max. power consumption: 30W
Maximum output current: 2.45A

POWER SUPPLY UNIT 11B

Model: TR30RDM150

Power supply (Volt): 15VDC
Operating temperature: 0°C – 40°C
Storage temperature: -40°C - +75°C
Max. power consumption: 30W
Maximum output current: 2.4A

13. LEGEND

01 Chest drain fixing
02 Keyboard
03 Display
04 Filter connection
05 Connection for intrapleural pressure sensing hose
06 Connection for power supply
07 Power ON/OFF
08 Suction line adjustment
09 Screen change
10 Liquid leak detection sensor
11A-B: Power supply unit
12 Power cord
13 Electrical plug (available in different versions)

Non-harmonised symbols used in Drentech®SilverLine:

SYMBOL	DESCRIPTION
	High suction pressure operation symbol (a)
	Charging in progress symbol Power supply connected with unit switched on (b)
	Battery status bar symbol (c)
	Charging in progress symbol / power supply connected with unit switched off (d)
	error symbol (e)
	Padlock symbol (keypad lock) (f)
	DO NOT DISPOSE (symbol g)
	over-tilt symbol (h)
	sound alarms disabled symbol (i)
	“patient” symbol (l)
	“reporting a problem” symbol (m)

Date of last version:
see last page: (REV: XX-XXXX)

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