
REDAX™

DRENTech™ SIMPLE

IT

ISTRUZIONI D'USO

DESCRIZIONE

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EN INSTRUCTION FOR USE

DESCRIPTION

The systems called Drentech® Simple are thoracic drains for adult/paediatric use fitted with a water seal valve separate from the collection chamber and designed for use by suction or gravity. Use by suction is possible by connecting the system to the portable vacuum units of the Drentech® line. The devices of this family are available in different versions based on the user's requirements.

The devices are composed of the following parts:

1. A collection chamber (1) made up of three graduated sections. The collection chamber has a maximum capacity of 2200 cc.
2. A water seal valve (2) that prevents air backflow to the patient's chest. The water seal valve has a graduated scale (3). The ball float valve (4) protects the water seal valve in the presence of high negative pressure preventing fluid backflow.
3. Some versions of the device are equipped with a dedicated valve located in the vacuum unit fastening connector (5). This valve has a dual function: it allows PALM EVO to correctly measure the intrapleural pressure in the receptacle and it acts as high-negative pressure relief valve when Simple is used by gravity.
4. The patient tube (6) is equipped with a lock-in connector (7) and a sliding clamp (8), which allow easily replacing the collection chamber, and an anti-kinking profile without a steel spring thus guaranteeing full compatibility with the magnetic fields generated

by Magnetic Resonance (NMR).

5. The handle with integrated hooks (9) facilitates transporting the drainage unit and also allows hooking the device onto the patient's bed.
6. Each system is equipped with a rear needle-free sampling point (12) for qualitative analysis of the drained fluid.

INSTRUCTIONS FOR USE

The Simple systems are suitable for: evacuating air and/or fluids from the thoracic cavity or the mediastinum; preventing build-up of air and/or fluids in the thoracic cavity; re-establishing and maintaining normal intrathoracic pressure gradients; facilitating complete pulmonary expansion to restore normal respiratory dynamics. The systems can perform their function by gravity or suction.

Note: This device can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

GENERAL WARNINGS

- The product is sterile if the package is intact. Do not use if the package has been opened or damaged.
- Dispose of the device after each single use; do not reuse.
- Avoid exposure to high temperatures and ultraviolet rays during storage.

WARNINGS/PRECAUTIONS

- The drainage system must be held in vertical position, including during transport, and always below the fluid level of the pleural cavity (at least 50 cm).
- Avoid looping or squashing the patient drain tube.
- Do not obstruct or cover the water seal chamber filling connector (10).
- Thanks to their particular design, the Simple systems are stable and do not require an IV pole to position them on the floor.
- Operation by suction is permitted only when connected to the Redax portable vacuum units which are expressly declared compatible. Suction can be applied by means of an external line provided that it is equipped with a thoracic regulator and under direct responsibility of the attending physician.
- The fluid contained in the collection chambers of the drainage system must be disposed of and may absolutely not be used for infusion.

PREPARING THE DRENTECH SIMPLE SYSTEM

The system is packed in double sterile wrapping. In order to use the product, take it out of its wrapping using a sterile technique.

a. Filling the water seal chamber

1. Prepare a preferably 60ml syringe.
2. Insert the syringe in the connector (10) and fill the seal valve (2) up to the broken line (about 45ml as also indicated on the front of the device). The water will turn blue thanks to a dye contained in the valve.
3. When the device drains by gravity, the water seal acts as pressure gauge. On the graduated scale of the relative column (3), you can qualitatively observe the intrathoracic pressure oscillations synchronous with the patient's respiratory acts.
4. Air bubbles may appear in the chamber underneath the water seal; these indicate air leakage from the pleural cavity, whose entity can be qualitatively measured (see "e. Patient air leakage diagnostics - point 1").

b. Connecting the patient tube

1. Close the sliding clamp (8) and connect the patient tube (6) to the thoracic catheter of the patient by means of the connector (11). Subsequently, connect the patient tube to the collection receptacle by means of the lock-in connector (7). Fit the connector until you hear a click indicating that it is correctly positioned. At this point, the system is ready to operate by gravity.
2. For operation by suction, prepare a Redax portable vacuum unit following the instructions given in the relative instruction leaflet. Connect the unit to the dedicated connectors and start using the device by controlled suction. Specific versions of the vacuum unit allow objectively measuring some clinical parameters, such as patient air leakage, etc.
3. To connect the device to the centralized suction, prepare a thoracic pressure regulator with an adjustment range within the safety parameters for the intended application. After adjusting the suction, connect it to the connector (10).
If you intend to use the system to treat neonatal patients, use the dedicated tubing for connection to the built-in Heimlich valve (code Redax 10564).

c. Disconnecting the patient tube

This operation can be carried out to replace the collection receptacle or to better arrange the patient tube, for example, based on the bed configuration.

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1. Close the patient tube by means of the sliding clamp (8).
- WARNING:** Never disconnect the quick-coupling before having closed the clamp.
2. Disconnect the connector (7) by pressing on the sides and at the same time pulling at the two ends. Then carry out the desired operations, for example, replace the collection receptacle.
3. To refit the connector, insert the two ends and press until you hear a click.

d. Resetting high negative pressure

During operation by gravity (this does not occur using the vacuum unit), high negative pressures may be created in the collection receptacle. These pressures spontaneously balance out over time. If you want to manually reset the pressure, operate as follows:

*** Systems equipped with a reset valve (5)**

In this case, just act on the dedicated valve using a ballpoint pen or similar. Pressing the valve will cause it to open and, consequently, the liquid in the graduated column of the water seal valve to go down.

WARNING: Stop the reset when the liquid level has reached zero. Do not keep the valve open for a prolonged time during operation by gravity because of the risk of pneumothorax.

*** Systems not equipped with a reset valve**

In this case, reset the pressure by disconnecting the lock-in connector of the patient tube (7) after closing the sliding clamp (8). Follow the procedure as described above under c.

e. Patient air leakage diagnostics

1. Air bubbles in the lower chamber of the water seal valve (2) are an indication of air leakage and, consequently, incomplete lung expansion. A reduction in the amount of bubbles over time indicates a gradual reduction of the leakage and better lung expansion. When lung expansion is complete, there will be no bubbles in the water seal valve.

WARNING: If the phenomenon persists, ask a physician to intervene after carefully checking the connections.

2. Some versions of the vacuum unit allow digital air leakage measurement (e.g. Palm EVO). Refer to the relative user manual.
- NOTE:** Evaporation, especially at high temperatures, or considerable air leakage from the patient, may cause the liquid level in the water seal valve to drop. It is advisable to periodically check the liquid level and top it up if necessary.

f. Transporting the device

The device should preferably be transported in vertical position and below the level of the patient. If necessary, the device can be placed in horizontal position as shown in Fig. 1. In this position, avoid mixing the liquids between the various chambers.

WARNING: In horizontal position, the water seal is compromised and you therefore need to clamp the patient tube and check that this condition is not excessively prolonged.

DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

CONTRAINDICATIONS

There are no known contraindications.

SINGLE-USE DEVICE. Reuse may lead to alterations in performance and cross-contamination risks.

STERILE - Sterilized by ethylene oxide MATERIALS USED

See the product technical data sheet.

MEDICAL DEVICE NOT MADE FROM NATURAL LATEX PVC-FREE PACKAGING.

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