



**REDAX™**

**EN**

## FLAT FLUTED DRAIN

### INSTRUCTIONS FOR USE

#### DESCRIPTION

The FLAT FLUTED drain, in pure biocompatible silicone, consists of:

- a flat part, in radiopaque white silicone, provided with 4 grooves which allow for effective drainage due to the capillarity principle;
- a part of transparent silicone tube, connected to the radiopaque end following an atraumatic profile, through which it is possible to connect to closed drainage systems by gravity or by suction.

Depth markings, placed at the end of the radiopaque end, which facilitate placement and removal of the drain. All drains are supplied in different sizes and are available with or without a steel trocar needle. All the configurations have the packaging in double sterile wrapping and contain a special fitting for the connection to the collection devices.

#### INTENDED FOR USE

The drains are intended for short-term use, in clinical cases, which require the evacuation of liquids accumulated in natural and/or newly formed cavities following surgery. They are suitable for use in conjunction with gravity or suction drainage systems. Depending on the size of the tube, this device can be used in adult, pediatric and neonatal patients.

**Note:** These devices can be used in safe conditions only in the applications and in the manner indicated in this instruction sheet, in accordance with the type of product itself. The manufacturer denies all liability arising from misuse and any use different from that indicated herein. **The procedure must be performed by trained personnel, specialized in anatomical landmarks, safe technique, and potential complications.**

#### SUGGESTIONS FOR USE

##### How to place the drain

1. Before placing the drain, carefully choose the type, location and size needed. These parameters depend on the type of surgery, the patient's physiology and the surgeon's experience.
2. Using aseptic technique, carefully remove the sterile chest drain from the protective wrap.
3. Place the drain into the wound following an aseptic procedure. In products equipped with a trocar needle, use the latter according to the usual positioning technique, proceeding from the inside of the wound outwards.
4. Use the depth markings on the drain for correct positioning.
5. Secure the drainage with a skin stitch or band-aid.
6. Connect the positioned tube to the selected drainage system according to the methods indicated in the respective instruction sheet. In the case of drains equipped with a Trocar needle, cut the Trocar needle and dispose of it before connecting to the drainage system. If necessary, it is possible to use the fitting included in the package, in the models where it is provided, to make the connection to the drainage system.
7. It is recommended to secure the connector connection to the drain tube and also to the catheter with tape.

##### Drain removal

1. Remove the tape or suture or any other type of fixing system.
2. Disconnect the drainage system.
3. Remove the drain by pulling it gently and immediately apply a dressing to the wound.

#### CONTRAINDICATIONS

There are no known contraindications.

#### GENERAL INSTRUCTIONS

The product is sterile if the packaging is intact. Single use, disposable product. Reuse can lead to performance alterations and risks of cross-contamination. Dispose of after each single use, do not reuse. Avoid exposure to high temperatures and ultraviolet light during storage. For the disposal and dismantling of the device, it is necessary to take appropriate precautions and comply with the legal

provisions in force on biological hazardous waste.

#### WARNINGS / PRECAUTIONS

1. Before inserting the drains, identify the insertion point and choose the correct size.
2. The choice of drain size is based on the type of pathology and surgery and on the experience of the surgeon.
3. It is absolutely forbidden to tamper with the drain and/or make additional holes with any tool. Only the reduction of the length of the fluted section is allowed, provided that the operation is carried out with a clean cut, without tears or burrs, and that the remaining section is not less than 15 cm, in order not to reduce the drainage capacity.
4. After removal, ensure that the drain is intact as this is an essential condition to prevent any drainage fragments generated by accidental lacerations from remaining inside.
5. It is recommended not to use the drain for a period longer than 29 days.
6. After positioning the device, ensure that the fluted portion is fully embedded within the wound and that the connections are tight.
7. The device must be removed by hand, pulling gently and avoiding sudden maneuvers. Avoid using metal tools such as clamps or devices that could break the device.
8. It is not recommended to use drains with high vacuum bottles.
9. If the device remains on site for very long periods, it may become difficult to remove. Be extremely careful during removal operations.
10. Do not use roller clamps or other metal devices, abrasive tools (such as cotton balls or gauze) for silicone tubes milking. The use of these tools can cause damage to the drainage surface by reducing its mechanical strength.
11. Do not use any alcoholic or aggressive liquid to lubricate the surface of the tubes before milking procedure.
12. Redax strongly recommend using Vaseline cream or oil or any other non-aggressive lubricant to facilitate tube removal.
13. Depending on their CH/Fr, catheters are compatible with various connectors available on the market (in addition to those present in the package). Be sure to check compatibility with the accessory you want to use before using it on a patient.
14. All connections must be firmly taped to minimize the risk of accidental detachment.
15. Properly dispose of sharps in sharps container in accordance with state/Osha standards for blood borne pathogens and/or institutional policy.

STERILE - Ethylene oxide sterilized

#### MATERIALS USED

Biocompatible and Haemocompatible silicone.

**MEDICAL DEVICE NOT MADE WITH NATURAL RUBBER LATEX**



Date of issue of the last version:  
see : (REV.: XX-XXXX)

**CE** 0123



Redax S.p.A.  
Via Galileo Galilei, 18  
Poggio Rusco (MN) Italy

REV.: 07-2021

COD. 00210