REDAX®

IT

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#### INSTRUCTIONS FOR USE

#### DESCRIPTION

P-Eggy Drain chest drain valve is a disposable device consisting of a medical grade silicone collection chamber with 100 ml collection chamber, needleless access port for fluid withdrawal from the collection chamber by means of syringe (on the bottom of the valve). Compression of the collection chamber makes it possible to apply moderate and temporary negative pressure onto the patient and check of the intensity of air leak.

The device consists of the following parts (Fig. 1):

- Medical grade silicone collection chamber with a connector (1) on top for direct connection of drainage tubes or by means of adapters (6).
- Inlet anti-reflux valve (2) to avoid potential reflux of drained fluids towards the patient.
- Silicone one-way seal valve (3) and air leak well (4) for the evacuation of the patient's air and its confirmation.
- Needleless sampling point (5) for the removal of liquid that may be present in the collection chamber by means of syringe.
- Adapters (6) to allow connection to set of various types and sizes.

#### CLINICAL CONTEXT

In cases of traumatic, spontaneous, iatrogenic pneumothorax, etc. it is necessary to apply a drainage that primarily allows the evacuation of air from the patient's chest cavity, preventing reflux. For this reason, a drainage catheter connected to a one-way valve (Heimlich or similar) and collection systems such as bags, drainage systems, etc.) is often applied to collect any drained liquid, usually a small amount. The patient with pneumothorax may be able to walk and in many cases is discharged at home with drainage until pneumothorax is reduced and the drainage tube is removed; for this reason, easy to mobilize and manage devices such as percutaneous drainage connected to the P-Eggy Drain thoracic drain valve are preferred.

#### INTENDED USE

The device is indicated for the evacuation of air from the chest cavity of the patient affected by pneumothorax of various origins (e.g. spontaneous, iatrogenic, traumatic, etc.) or in cases in which patient mobilization with drainage is required. This can restore normal intrathoracic pressure gradients and facilitate complete lung expansion, restoring normal respiratory dynamics. This device can be

used in adult and pediatric patients.

#### GENERAL WARNINGS

- The product is sterile if the packaging is intact. Do not use if the package has been opened or damaged.
- Avoid exposure to high temperatures and ultraviolet rays during storage.
- Dispose of after each single use, do not reuse or re-sterilize.
- Do not clog the seal valve (03).
- Do not clamp the patient's tube during use.
- Do not use or prick the needleless access port with a needle.
- Do not leave the syringe in the needleless sampling port.

**NOTE:** This device can be used in safety conditions only for the applications and with the methods indicated in this operating instruction sheet, according to the type of the product itself. The manufacturer declines all liability arising from misuse and any use different from that indicated. **The procedure should be carried out by trained and experienced personnel and aware of safe techniques and potential complications.**

#### WARNINGS/PRECAUTIONS

- **The device should be positioned and used by specialized personnel, aware of the implications of the chest drainage procedure.**
- **The device shall be handled by wearing protective gloves and all devices specified in hospital protocols to prevent cross-contamination of the patient and/or operator.**
- **Any liquid contained in the collection chamber must never be reinfused, but disposed of with the drainage system in accordance with the current rules on biologically hazardous waste.**
- **Replace or empty the P-Eggy Drain chest drain valve when the volume of collected liquid exceeds the maximum capacity.**
- **If damaged, the P-Eggy Drain chest drain valve should be replaced.**
- **If there is evidence of occlusion, replace the system.**
- **Do not make a Y-connection between two catheters and a P-Eggy Drain valve.**
- **Make sure that the P-Eggy Drain valve is securely connected to the catheter.**
- **Make sure the system is upright to ensure proper drainage of patient's fluids.**
- **Make sure the device does not remain squeezed/pressed after emptying the collection chamber to ensure proper drainage of patient's fluids.**
- **Make sure the silicone collection chamber is free from fluids before compressing the device; otherwise, proceed with its emptying starting from the needleless sampling point, using a syringe**
- **The connectors, present in the package or the kit of compatible fittings (code 10543), allow connecting different types of drains to the P-Eggy Drain chest drain valve, as shown in Fig. 2 and Fig.3. Be sure to check the compatibility with the catheter you want to use before you use it on the patient.**

#### PREPARATION FOR USE

Proceed as follows:

1. Pull the device out of its sterile casing.
2. Visually verify that the device is intact and looks like in Fig. 1 and that all accessories provided are present.
3. Connect the drainage to the conical connector (1) of the P-Eggy Drain valve using the adapters in the package as shown in Fig. 2. Make sure that the connection is firmly secured.
4. Place the system, prepared as described, below the wound level to immediately begin drainage by gravity.


#### Patient drainage

The device allows drainage by gravity of air and liquids from the chest cavity of the patient. Slow and gradual compression of the silicone collection chamber allows applying a moderate and temporary suction (from -25 cmH<sub>2</sub>O to -27 cmH<sub>2</sub>O) and checking that the drainage tube is opened.

#### Sampling and emptying of the device

If the collection chamber is almost full, the liquid can be removed by connecting a syringe, preferably equipped with a luer-lock connector

Fig.3b

CH-Fr 26-28

Fig.2



Fig.3a

CH-Fr 20-24

to the needleless sampling point (5). This procedure should be carried out before compressing the collection chamber, according to the current hospital provisions on infection prevention. Clean with alcohol or disinfectant solution available the needleless port before and after applying the syringe.

**Note:** accidental leak of fluid due to system overturn or squeezing is usually prevented by the safety valve installed on the system.

**Patient air leak confirmation**

It is possible to confirm the patient's air leaks by adding about 1 ml of water into the air leakage well (4). The appearance of bubbling in water (air bubbles) is an indicator of air leak from the patient. Reduction in the amount of such bubbles over time indicates a progressive reduction in leak and lung expansion. Empty the air leak well after use. Alternatively, it is enough to slowly compress the silicone collection chamber and visually check the expansion of the same: the faster the expansion will be, the greater the air leak. When the lung expansion is complete, the collection chamber will remain compressed.

**Disconnecting the system**

Clamp all catheters before disconnecting the chest drainage valve from the patient.

**REMOVAL/DISPOSAL**

For the disposal and dismantling of the device, it is necessary to take appropriate precautions and comply with the legal provisions in force on biologically hazardous waste.

**CONTRAINDICATIONS**

There are no known contraindications.

**DISPOSABLE PRODUCT.**

Reuse can lead to performance alterations and risks of cross-contamination. Dispose of after each use.

Do not reuse or re-sterilize.

**STERILE - Ethylene oxide sterilized**

**MATERIALS USED :** Silicone, PP, ABS, PVC.

**MEDICAL DEVICE MADE WITHOUT NATURAL LATEX**

**Reference**

- 01. Connector for drainage tube
- 02. Inlet anti-reflux valve
- 03. One-way silicone seal valve
- 04. Air leak well
- 05. Needleless sampling port
- 06. Adapters

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