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#### 1. DESCRIPTION OF THE DEVICE

ALIVIA is a percutaneous drainage device available in a range of sizes. The product consists of the following components:

1. A main body with three inlets (Fig.1).
2. The first inlet (1) on the back is reserved for the introduction of the Veress needle. The two remaining inlets (2) are equipped with luer-lock connections and allow the connection of drainage devices such as syringes and/or collection systems, as described in this instruction sheet. On the main body there is a selector switch (3) which allows, when rotated by 90°, for manual drainage (pull&push) by gravity or suction. The drainage catheter is attached to the main body (4) and features holes and depth markings.
3. There is a Veress needle (5) for introducing the catheter. This needle has a mandrel with an atraumatic tip that protrudes from the end of the needle. This mandrel can slide longitudinally until the cutting profile of the needle itself is uncovered during insertion, and then returns to its extended position by means of a spring. The needle is equipped with a handle (6) in which there is a visual indicator that informs the user of the position of the protective atraumatic mandrel.
4. The catheter is made of biocompatible, pyrogen-free and totally radiopaque material, allowing easy localisation by X-ray.

#### STERILE - STERILIZED WITH ETHYLENE OXIDE

The product is sterile if the packaging is intact. The device is single use and therefore requires no cleaning or maintenance.

#### STORAGE

The device must be stored according to the label: temperature limit (Min/Max) 0°C/60°C.

#### 2. DIRECTIONS FOR USE

The ALIVIA system is primarily intended to provide access, in a minimally invasive manner, to the thoracic and abdominal cavities of the patient in order to evacuate fluids or air accumulated there, or a combination of the two. The device can be used for gravity or suction drainage, connected to collection systems.

#### PATIENT GROUP

The device is indicated for use on adult patients.

#### USER GROUP

All products are intended for use by qualified medical personnel, such as doctors and nurses.

#### MEDICAL CONDITIONS

The areas of application of ALIVIA drainage systems are treatment of pneumothorax (primary and secondary), malignant or chronic pleural effusions, empyema and abdominal ascites.

#### CONTRAINDICATIONS

Coagulopathies and other blood diseases

#### WARNINGS

- Do not re-sterilise the device.
- Single Use. Reuse can lead to altered performance and risk of cross-contamination.
- Do not use the device beyond the expiry date indicated on the packaging.
- Avoid exposure to high temperatures and direct UV light during storage.
- Tampering with the device and/or its parts is forbidden. This device may only be used safely in the field of application and in the manner specified in this instruction sheet, corresponding to the product type. The manufacturer denies all liability arising from misuse and any use different from that indicated herein. Check that the visual indicator of the mandrel position runs freely and that the corresponding indicator is green when the mandrel is in the extended position. If, during initial compression, the green indicator on the outer surface of the chest/abdominal wall does

EN

## MULTIPURPOSE DRAINAGE CATHETER

INSTRUCTIONS FOR USE

#### DECLARATION OF LIABILITY

This document is the exclusive property of Redax S.p.A. all rights reserved. Any disclosure, reproduction or transfer of content to third parties is prohibited without the prior authorization of the Company. Redax S.p.A. denies all liability with respect to damage to people or property due to improper use of this product and failure to observe the indications, warnings, instructions and precautions in this user manual. In the event of a serious accident related to the medical device, the user must notify the manufacturer and the competent authority of the state in which the event occurred.

This user manual is supplied in printed form and accompanies the device.

not disappear, do not continue needle introduction and replace the device.

- The device is designed to be used with Redax products and accessories. If alternative products and accessories are used, it is the doctor's responsibility to check their compatibility before use.
- Do not use the device for more than 30 days.
- The duration of the drainage procedure must be determined by the doctor based on clinical considerations.
- Do not use the device if the casing or any of its parts are damaged.
- Do not manipulate the drainage in any way or drill additional holes.
- In the case of suture fixation, carefully avoid suturing the tube itself or cutting it with sharp needles. Avoid excessively tight sutures as this may cause drain rupture.
- After positioning the drain, ensure that the perforated section is completely included within the wound and that the connections are tight.
- Drain should be removed manually by gently pulling, avoiding sudden manoeuvres.
- Avoid using metal tools such as clamps or devices that could break the drain.
- After removal, check that the drainage set is complete as this is essential in order to exclude the presence in the cavity of drainage fragments caused by accidental tearing.
- Improper use could expose the patient and the operator to various dangers.
- Keeping the drain *in situ* for very long periods could make its removal difficult.
- Care should be taken during removal operations.

### 3. CLINICAL BENEFIT

The clinical benefit associated with ALIVIA percutaneous drains is to create a minimally invasive access to the patient's thoracic or abdominal cavities in order to evacuate accumulated fluids (air and pathological fluids) and alleviate related symptoms. The device's main feature consists of an integrated dual one-way valve that ensures patient protection from fluid ingress and simplified user control of the device, which allows for: gravity drainage and controlled manual suction (pull & push manoeuvre), all without the use of additional components resulting in improved device management.

### 4. INSTRUCTIONS FOR USE

#### PREPARATION FOR USE

Before using the system, carry out the preparation procedure as described below (Fig.2):

1. Check the intact state of the packaging to ensure the functionality and sterility of the product.
2. Open the package using sterile procedure
3. Determine the position for drainage insertion with the aid of an ultrasound or radiological investigation.
4. Prepare a sterile field around the chosen area, taking normal precautions for skin disinfection.
5. Remove the components required for use from the packaging, such as the main body of the device with the attached catheter and Veress needle.
6. Insert the needle through the rear access (1) until the handle is fully retracted. The needle and its protective mandrel should protrude a few millimetres from the end of the catheter.
7. Perform local anaesthesia on the site selected for needle introduction using the usual anaesthetics.
8. If the doctor considers it necessary, in order to facilitate the introduction of the catheter, make a small incision in the skin. The scalpel included in the package can be used for this purpose.

#### Valve hub management

**Initial position (longitudinal):** In this position, the needle for drainage placement is inserted.  
The drainage is closed.

**90° rotation:** In this position, it is possible to drain fluids by gravity, suction or manual technique (push & pull).  
The drain is open, but the patient remains protected by

the internal one-way valve, which prevents any reflux.

#### NEEDLE INTRODUCTION AND CATHETER PLACEMENT

After carrying out the preparation described above,

1. Hold the device firmly, using both hands. The hand located on the main body is responsible for providing the insertion thrust, directing the needle and keeping the needle handle in contact with the body of the device. The other hand is recommended to be placed close to the distal end; grasping the catheter with the thumb and forefinger will conveniently avoid inserting the needle too far into the patient's pleural cavity once it has passed the skin.
2. Place the needle tip against the patient's skin and apply gentle pressure until the green visual indicator disappears (Fig.3).
3. Introduce it very slowly until the green visual indicator reappears; this means that the end of the needle has entered the pleural space and the protective mandrel has returned to its extended position, fulfilling its safety function (Fig. 4).
4. Tilt the handle of the device so as to orientate the direction of catheter introduction. To verify that the introduction has taken place at the desired point, i.e. at the area where the fluid to be drained is accumulated, exploratory sampling can be carried out. To do this, connect the syringe to the luer-lock connector on the needle by removing the corresponding cap, retract the plunger and visually observe (Fig.5) the fluid entering the syringe.
5. Begin removal of the needle and, at the same time, complete catheter placement, as shown in Fig.6.

#### FLUID EXTRACTION: (PULL&PUSH Fig.7)

Once the catheter has been positioned, draw out the fluid to be drained with the aid of a syringe(7) and collection bag (8). To connect these accessories, open the corresponding plugs (2) located on the two luer-lock connections on the main body. Connect the syringe to one of the two connectors and the collection bag to the other.

The choice of which connector to use to connect the syringe or bag is left up to the doctor, who will decide on the basis of ease of operation. It is necessary to turn the selector switch on the main body clockwise as far as it will go, corresponding to a 90 degree angle to the initial position. Next, pull the syringe plunger (Fig.7a): the fluid will be sucked from the chest cavity into the syringe. By pressing the plunger, the fluid will automatically be sent to the collection bag (Fig.7b). Repeating the pulling and pressing operation in sequence will result in the evacuation of the fluid, which will be collected in the bag. At any time, the bag can be replaced or the syringe or both can be removed, depending on the outcome of the drainage. To do this, it is first necessary to reposition the dial longitudinally (knob oriented according to the axis of the catheter). It is advisable, when removing the accessories, to close the protective caps to avoid possible contamination. In any case, if the caps remain open, there is no risk of reflux towards the patient.

#### GRAVITY DRAINAGE: COLLECTION BAG

After the catheter has been positioned and secured, it is possible to proceed with the drainage of fluids by gravity. In this case, the drainage bag included in the package, through which drained liquids can be collected by gravity, must be connected to any of the available fittings (02). Correct and safe operation of the device is only guaranteed if Redax-compatible devices are used, as indicated in Section 7 of the current operating instructions. The use of non-compatible components may compromise the performance and safety of the system.

To drain the fluids, proceed as follows (Fig.8):

1. Connect the collection bag to one of the luer-lock fittings (02). Before carrying out subsequent operations, make sure that the other fitting is closed by the appropriate plug.
2. Place the collection system below the level of the patient.
3. Turn the selector switch clockwise as far as it will go (Fig.8a). At this point, there will be a connection between the catheter and the collection system with an interposed non-return valve.
4. If you wish to replace the collection bag, e.g. in case of complete filling, turn the selector switch to the home position and then replace the bag. After replacement, turn the selector switch 90° again.

#### FLUID DRAINAGE: POST-SURGERY COLLECTION SYSTEM

After the catheter has been positioned and secured, it is possible to proceed with suction or drop fluid collection using an appropriate postoperative collection system (e.g. silicone reservoir, chest

FIG. 9a

FIG. 8

**GRAVITY DRAINAGE**



FIG. 8a

drainage system, etc.).

To connect a postoperative drainage device, proceed as follows:

1. If a chest drainage system is to be used, connect the connector of the collection system used to one of the luer-lock fittings (2).
2. Prepare the collection system by carrying out the operations prescribed by the manufacturer, for this purpose refer to the operating instructions of the respective products; then connect the system to the free end of the connector. Full compatibility with Redax drainage systems such as Drentech WearABLE (Fig.9a) and P-Eggy drain (Fig. 9b). In this case, use the fittings included in the respective packages.
3. Turn the selector switch, located in the main body, clockwise as far as it will go, corresponding to a 90° angle to the initial position. At this point, there will be a connection between the catheter and the collection system with an interposed non-return valve.
4. If you wish to replace the collection system, e.g. in case of complete filling, turn the selector switch to the home position and proceed with the replacement. When finished, turn the selector switch 90° back to the previous position.

**CATHETER FIXING**

A fixing system is interposed between the body of the ALIVIA system and the catheter to ensure stable securing to the patient's skin by means of sutures (10).

**DEVICE REMOVAL**

At the end of the drainage operation, the device can be removed. For removal, proceed as follows:

1. Remove any bandages and dressings and disinfect the area around the point of entry of the drainage tube.
2. Remove the fixing suture and prepare the wound closure suture (e.g. tobacco bag) in advance.
3. Pull the drain tube out by exerting moderate traction. At the end of the removal, check the intact state of the tube and then proceed to close the wound.

**5. PACKAGING CONTENT**

The packaging consists of the following components:

- polyurethane catheter with integrated valve system (Fig1)
- Verres needle (6)
- 2000 ml collection bag (8)
- one syringe (7);
- one disposable scalpel (9);

**6. COMPATIBLE ACCESSORIES**

The compatible accessories made by Redax are:

- Luer-lock connector code 10504
- 2000 ml drain bag

For codes, please contact [info@redax.it](mailto:info@redax.it).

**7. COMPATIBLE DEVICES**

The ALIVIA system is compatible with the following devices manufactured by Redax:

- Drentech WearABLE
- P-Eggy drain: chest drain valve

For codes, please contact [info@redax.it](mailto:info@redax.it).

**8. DISPOSAL**

When disposing of the device after use, appropriate precautions must be taken and the current legal provisions concerning biologically hazardous waste must be observed.

**9. MATERIALS USED**

- Polyurethane (PU) for the drain part in contact with the patient.
- PVC integrated fixing system

**LATEX- AND PHTHALATE-FREE MEDICAL DEVICE**

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