
FIG. 3b

FIG. 4

CE 0123

Redax S.p.A.
Via Galileo Galilei, 18
Poggio Rusco (MN) Italy

REDAX™

IT

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INSTRUCTIONS FOR USE

DECLARATION OF LIABILITY

This document is the exclusive property of Redax S.p.A. all rights reserved. Any disclosure, reproduction or transfer of content to third parties is prohibited without the prior authorization of the Company. Redax S.p.A. denies all liability with respect to damage to people or property due to improper use of this product and failure to observe the indications, warnings, instructions and precautions in this user manual. In the event of a serious accident related to the medical device, the user must notify the manufacturer and the competent authority of the state in which the event occurred. This user manual is supplied in printed form and accompanies the device.

MANUFACTURER


Redax S.p.A.
Via Galileo Galilei, 18, 46025 Poggio Rusco (MN) Italia
VAT NO: 02556750368
info@redax.it
www.redax.it

1. DESCRIPTION OF THE DEVICE

The Drentech WearABLE chest drain is a single-use device consisting of an air-liquid separation chamber equipped with an anti-reflux valve and two one-way valves for air evacuation connected to a 300 ml collection chamber.

The device consists of the following parts (Fig.1):

- Air-liquid separation chamber (01) fitted with a multi-calibre connector (02) on top for connecting drainage pipes. This chamber is equipped with a hydrophobic filter, to ensure the separation between air and liquids.
- One-way valve located at the inlet of the separation chamber (03) to prevent potential reflux of fluid to the patient.
- One-way silicone low-pressure relief valve (04), positioned on the air-liquid separation chamber, which allows the immediate expulsion of any air from the patient's chest cavity, preventing it from re-entering. The valve is equipped with a special reservoir (pocket) (05) which, when filled with a few drops of saline or water, provides a visual indication of the patient's air leakage by bubbling.

EN

- Silicone safety valve for evacuating air in the event of high flow rates (06), such as those resulting from patient coughing or temporary unavailability of the primary valve (04).
- Button to facilitate the transfer of liquid from the collection chamber to the bag (07). This button also acts as a negative pressure release valve.
- Flow indicator (15) to perform an evaluation of the negative pressure in the system, to check the patency of the drain and also to assess the presence of air leak from the patient.
- Flexible collection reservoir with a capacity of approx. 300ml (08) and equipped with a needle-less valve (09) to enable it to be emptied in order to restore the original collection capacity. On the sides, there are two slots (10) to allow the device to be attached to an anatomical belt (11) in the manner described below (Fig.3a), allowing full mobility. Alternatively, a handle (12) is included for securing the device to the patient bed, as described in Fig.3b.
- Rubber patient connecting tube (13), 500 mm long, fitted with a CH 28-36 connector. Three further connectors are also included in the package: a 15 to 22 CH multi-calibre connector (13a), a 24 to 28 CH multi-calibre connector (13b) and a luer-lock cone for connecting the WearABLE to Drentech Unico (13c - 13d) percutaneous drainage catheters.
- 600 ml collection bag (14) with liquid drain stopcock and male luer-lock connection for connection to the collection container. The bag is equipped with an integrated "flutter" valve to prevent the return of liquid to the container and a 500 mm connection line.

STERILE - STERILIZED WITH ETHYLENE OXIDE

The product is sterile if the packaging is intact.

The device is disposable and therefore requires no cleaning or maintenance.

STORAGE

The device must be stored according to the label: temperature limit (Min/Max) 0°C/60°C.

2. DIRECTIONS FOR USE

The device is indicated for the collection of fluids such as air, liquid or a mixture of both accumulated in the thoracic cavity

PATIENT GROUP

The device is indicated for use on adult and paediatric patients.

USER GROUP

All products are intended for use by qualified medical personnel, such as surgeons and nurses.

The operations listed in the paragraph "instructions in an out-of-hospital environment" are aimed at non-expert users, such as patients (only if adults) or any caregivers, suitably trained by qualified medical personnel.

MEDICAL CONDITIONS

The fields of application for Drentech WearABLE chest drainage are: post-operative procedures following surgery and/or trauma, pneumothorax of various origins (e.g. spontaneous, iatrogenic, traumatic, etc.), recurrent and/or chronic pleural effusions.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS/PRECAUTIONS

- **Do not re-sterilise the device.**
- **Disposable product. Reuse can lead to altered performance and risk of cross-contamination.**
- **Do not use the device beyond the expiry date indicated on the packaging.**
- **Avoid exposure to high temperatures and direct UV light during storage.**
- **Tampering with the device and/or its parts is forbidden.**
- **This device may only be used safely within the intended area of use and in the manner specified in this instruction sheet. The manufacturer denies all liability arising from misuse and any use different from that indicated herein.**
- **The device is designed to be used with Redax products and accessories. If different products and accessories are used, it is the doctor's responsibility to check their compatibility before use.**

- **Do not use the device for more than 30 days.**
- **The duration of the drainage procedure must be determined by the doctor based on clinical considerations.**
- **Do not use the device if the casing or any of its parts are damaged.**
- **Improper use could expose the patient and the operator to various dangers.**
- **The liquid contained in the collection tank must be disposed of and under no circumstances should it be used for infusion.**
- **Out-of-hospital environment:**
- **In case of impact, accidental fall or tugging of the device, make sure it is intact and, if in doubt, contact the qualified reference staff.**
- **If the device gets disconnected, you cannot transfer in the bag and/or if there are any leaks, immediately contact the reference health facility and follow the instructions provided during discharge.**
- **In case of shortness of breath, fever, redness or localized pain at the drain insertion point or any other symptom indicated as relevant and significant by the medical personnel at the time of discharge, contact the reference health facility.**
- **If the drainage tube moves from the insertion point, immediately contact an ambulance service reporting the event and informing the operator of a potential risk of "pneumothorax".**

3. INSTRUCTIONS FOR USE

3a. INSTRUCTIONS FOR IN-HOSPITAL ENVIRONMENTS PREPARATION FOR USE

Before using the set, carry out the preparation procedure as described below:

1. Remove the device from its sterile casing.
2. Visually check that the device is intact, as shown in Fig.1 and that all required accessories are present.
3. Connect the chest drain coming from the patient to the CH 28-36 multi-calibre connector pre-connected to the 500 mm long rubber hose in the package (13) and then connect the hose to the Drentech WearABLE multi-calibre connector (02). If the patient's chest drain is smaller (15-22 CH or 24-28 CH), or has a luer-lock connection, use the fittings included in the package (13a, 13b, 13c and 13d). Make sure the connection is secure.
4. Place the system, prepared as described, below the level of the wound to immediately begin gravity drainage. The strap in the package can be used to maintain the upright position.

PATIENT DRAINAGE

The device allows gravity drainage of air and liquids from the patient's chest cavity. By means of the air-liquid separation chamber and its hydrophobic membrane (01), located immediately next to the catheter connection, it is possible to separate the liquid portion of the drainage from the air portion: the former is conveyed by gravity into the collection tank (08); the air is immediately expelled by the special one-way silicone valve (04).

CHECKING FOR AIR LEAKS FROM THE PATIENT


Patient air leaks can be confirmed by adding about 1 ml of water to the air leak reservoir ("pocket") (05) located on the one-way valve. The appearance of bubbles in the reservoir (air bubbles) is an indication of air leaking from the patient. The reduction in the quantity of these bubbles over time indicates a progressive reduction in leak and subsequent lung expansion.

The verification of air leaks can also be carried out by observing the re-expansion of the flow indicator (15). After compression of the collection reservoir, as shown in Fig.2, the indicator will collapse on itself; in the presence of a patient air leak, its progressive re-expansion will be observed over time. In particular, complete re-expansion is achieved with the entry of 60 cc of air into the container. The ratio between this volume and the time taken for re-expansion can provide a measure of the patient's air leak.

AIR LEAK EVACUATION : HIGH FLOW RATE

In the event of a sudden evacuation of air with high flow rates, e.g. during a coughing fit, the safety valve (06) enables rapid evacuation without dangerous overpressure. The simultaneous action of the two valves, (04) and (06), allows high flows to be drained in short

FIG. 1

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moments. The additional valve also has a safety function, preventing the drainage from blocking even if the primary valve is unavailable.

NEGATIVE PRESSURE APPLICATION

It is possible, at any time, to apply temporary, moderate suction by simply compressing the flexible tank (08), as shown in Fig.2.

Compression can generate suction down to -60cmH2O, the duration of which depends on the patient's leak flow. This may be performed to assess the patency of the drainage catheter or to facilitate drainage.

The flow indicator (15) allows evaluation of the presence or absence of negative pressure within the system. Following the aforementioned compression of the collection tank, the indicator will collapse on itself, indicating the presence of negative pressure, which will be gradually maintained until the complete re-expansion of the indicator itself.

FIXING THE DEVICE

The WearABLE drain is primarily designed to enable early mobilisation of the patient. For this, a 180 cm long anatomical belt (11) is available, which can be fastened as shown in fig.3a.

The length is adjustable and a additional belt is also available to purchase that can be fastened on the patient's hips or worn as a shoulder strap. Alternatively, during the first hours of the patient's hospitalisation, a fixing handle (12) is provided with which WearABLE can be attached, for example, to the patient's bed (fig.3b).

EMPTYING THE COLLECTION TANK

The collection tank (08) has a maximum capacity of about 300ml. At any time, this can be emptied in order to restore full collection capacity. To perform this operation, simply connect the appropriate bag to the needle-less valve (09), located at the front bottom of the reservoir itself. The liquid will transfer from the reservoir to the bag by gravity, as shown in Fig. 4

NOTE: transfer to the collection bag generally takes place spontaneously without the need for any intervention by the operator. In rare cases, in the total absence of air leakage and/or in the presence of negativity in the reservoir, an excessively slow transfer may occur. This condition is highlighted by the flow indicator, which appears to be partially or totally collapsed. In this case, pressing the button (07) allows a small amount of air to enter from the outside, which allows the liquid to flow down to the bag quickly. The presence of the filter membrane prevents direct contact of the tank with the external environment. After the transfer, the bag must be disconnected and disposed of. Make sure the tap is closed, use the white clamp to isolate the contents of the bag and proceed with disposal in accordance with the provisions in force on biologically hazardous waste.

SYSTEM REMOVAL

Clamp all catheters before disconnecting the Drentech WearABLE chest drain from the patient.

If the patient continues to be treated with Drentech WearABLE in an out-of-hospital environment, it is the responsibility of the reference medical personnel to provide all the indications on the device management, on how to behave and who to contact in case of symptoms potentially related to the drainage state or the device operation.

3b. OUT-OF-HOSPITAL INSTRUCTIONS

All operations relating to the positioning and start-up of the drainage procedure must be carried out in a hospital environment by expert and qualified staff.

DEVICE MANAGEMENT

- Keep the device below the insertion point of the drain.
- Check that the patient connection tube (13), if any, is well extended, not twisted and well connected.
- If the Drentech WearABLE device gets disconnected, reconnect it and contact the hospital or the contact details provided at the time of discharge.
- If it is necessary to replace, reposition or adjust the fastening belt, refer to the "fastening the device" paragraph of section 3a. If necessary, the hospital facility will be able to provide a new fastening belt.

PERSONAL HYGIENE

The device and the insertion point of the drain should not be

immersed in water, such as in a bath. It is preferable to wash the various parts of the body separately or take a shower. If possible, use the support of a caregiver (assistant) who can keep the device dry or perform any procedures indicated at discharge.

EMPTYING THE COLLECTION TANK

Follow the instructions recommended by qualified medical personnel on whether and when to empty the system in a bag and how to record the amount of drained liquid.

To empty the Drentech WearABLE collection system, after ensuring that the tap of the appropriate bag is closed, simply connect the bag to the needle-less valve (09), located at the front bottom of the reservoir itself. The liquid will transfer from the reservoir to the bag by gravity, Fig. 4

If the liquid is not transferred into the bag, the button (07) must be pressed, which allows a rapid emptying into the bag.

After the transfer, the bag must be disconnected and disposed of. To this end, isolate the contents of the bag using the integrated white clamp, disconnect the collection bag and, finally, proceed with disposal in accordance with the methods agreed with the reference healthcare facility and according to the regulations in force regarding biologically hazardous waste.

4. ACCESSORIES CONTAINED IN THE SET

The sets usually contain:

- Drentech WearABLE chest drainage device
- 1 x Strap to help mobilise the patient and maintain an upright position
- 3 x fittings
- Collection bag
- Patient tube
- 1 Fixing handle
- Leaflet

5. COMPATIBLE ACCESSORIES

- Connector supplied by Redax
- Belt provided by Redax
- Bag provided by Redax

For the accessory codes that Redax guarantees compatibility, please contact info@redax.it.

6. DISPOSAL

When disposing of the device after use, appropriate precautions must be taken and the current legal provisions concerning biologically hazardous waste must be observed.

7. MATERIALS USED

Cover: SBC

Air-liquid separator filter: hydrophobic acrylic copolymer membrane

Valves: silicone and synthetic rubber.

Liquid collection tank: PS

Hose for connection: rubber

Fittings: PP

Needle-less valve: PC

Leak connector cover: ABS

Handle: PP

Flow indicator: rubber

8. LABELLING AND SYMBOLS USED

In accordance with standard ISO 15223-1.

KEY

- 01 Air-liquid separation chamber
- 02 Multi-calibre connector
- 03 Anti-reflux one-way valve
- 04 Primary one-way valve for air evacuation
- 05 Air leak detection pocket
- 06 Additional high-flow air evacuation valve
- 07 Liquid transfer button and negative release valve
- 08 Flexible drained liquid collection tank
- 09 Needle-less valve for tank emptying
- 10 Loops for belt connection
- 11 Anatomical patient fastening strap
- 12 Fixing handle for bed/stand
- 13 Patient connection tube with 28-36 CH connector
- 13abcd: Connectors

- 14 Collection bag for drained fluid
15 Flow indicator

Date of issue of the latest version:
see last page: (REV.: XX-XXXX)

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