EMPTYING THE COLLECTION BAG

At the end of the fluid transfer into the bag, it must be disconnected and emptied.

How to do it :

1. Disconnect the bag by turning the connection counterclockwise, as shown in figure 1. The connection port is self-sealing, so when the bag is removed it prevents any fluid leaking from the semi-rigid container.
2. Stand over a sink. Open the blue tap at the bottom of the bag by sliding it in the direction shown by the arrow in figure 2 (towards the graduated scale side).
3. Once empty, close the bag again by sliding the blue tap in the opposite direction, so the bag is sealed and ready to be used again.
Note: please, follow the instructions on fluid disposal given by your doctor at discharge.
4. Finally, clean the system and all accessories using a neutral solution or detergent.

WHAT MAY HAPPEN

Accidental removal of the drain :

- **Do not** try in any way to put the drain back in place.
- Disinfect and cover the skin opening with sterile gauze and a post-operative dressing.
- Contact your doctor or the ward using the number shown on the first page of this guide

Fever ($\geq 38^{\circ}\text{C}$)

or redness around the drain :

contact your doctor or the ward using the number shown on the first page of this guide.

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WHAT CAN BE MONITORED

FLUID DRAINAGE :

To measure fluid drainage, use the graduated scale on the semi-rigid container. Check and record the amount of fluid collected as often as instructed by your doctor at discharge.

AIR LEAK :

Your doctor may ask you to check for air leaks. To do this there are two modalities:

- Add few drops of water into the front "chamber": bubbling means that an air leak is present.

- Create negative pressure in the system by compressing the semi-rigid container, following the instructions shown on the front panel. The monitor the yellow cylindrical indicator: depending on the time needed for the indicator to fully re-expand, the doctor will be able to assess the patient's air leaks.

TRANSFERRING FLUID INTO THE BAG

When to do it:

- Whenever the rigid container reaches 300 ml (full)
- Before going to bed
- As often as instructed by your doctor at discharge

Steps to follow :

1. Wash your hands.
2. Connect the collection bag to the special port at the bottom of the rigid container by pushing and turning clockwise at the same time, as shown in figure 2.
3. Wait for the fluid in the semi-rigid container to flow into the bag by itself. If this does not happen, press and hold the black button shown in figure 3 until all the fluid has been transferred.
4. Wash your hands.

PRECAUTIONS TO FOLLOW WHILE USING THE DEVICE

1. Always keep the system below the level where the drainage tube enters the chest. Fix it using the belt around your waist, wear it over the shoulder, or place it as shown in the picture below.
2. Avoid sudden movements that could cause the drainage tube to move out of place.
3. During rest or sleep, place the container next to the bed, always keeping it below the level of the chest tube insertion point. It is recommended to empty the semi-rigid container (using the steps described in the section "Transferring fluid into the bag") before going to bed.

AROUND WAIST

OVER SHOULDER

TWO BELTS TOGETHER
(provided by the doctor/available separately)

PERSONAL HYGIENE

The device is designed to be used even while taking a shower. Therefore, it does not need to be removed before.

After washing, simply dry the cotton belt or replace it with a clean, dry one.

When possible, it is recommended to have a family member or caregiver help you during the shower.

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To be completed by hospital staff before patient discharge

EN

NAME AND SURNAME:

HOSPITAL NAME :

REFERRING DOCTOR AND EMERGENCY CONTACT DETAILS :

Scan QR Code
to download
this document

Dear Patient,

this document explains what the chest drain used during your recent hospital stay is, what it is for, and how to take care of it at home.

WHAT IS IT FOR ?

Your lungs need enough space to expand when you breathe. If air and/or liquid collect in the space between the lung and the chest wall, the lung cannot work properly.

The device you received is a drainage system connected to a small tube placed in the area where the fluids collected. It helps remove the air and/or liquids from the chest and prevents them from going back inside. This allows the lung to return to normal function.

WHAT DOES IT LOOK LIKE?

BELT

COLLECTION BAG