




REDAX™

MEDICAL DEVICES



CONNECTORS

Connectors from 20 CH to 28 CH
and luer lock adapter

P-EGGY DRAIN Chest drain valve

Included components:
Chest drain valve and connectors

Code 10808

Packaging: 10 units per box

CONNECTORS KIT

Code 10543

Packaging: 10 units per box



Scan QR Code
to access
product video




P-EGGY
DRAIN
CHEST DRAIN VALVE



Via Galileo Galilei, 18 - 46025 Poggio Rusco (MN) ITALY - Tel. +39 0386 830582

info@redax.it | www.redax.it

REV. 06-2026



**SMALL and
PORTABLE
DRAIN VALVE**



P-EGGY
DRAIN
CHEST DRAIN VALVE



A LIGHTWEIGHT SOLUTION FOR
EASY PATIENT MOBILITY



ACTIVE DRAINAGE
Squeeze the reservoir to create suction.



VALVED SYSTEM
Double integrated one-way valve provides excellent protection.



AIR-LEAKS DETECTION
Thanks to the integrated well, air leakage may be detected by bubbling.



NEEDLELESS LUER LOCK CONNECTOR
port for fluid sampling and drainage discharge.

HOME-CARE



An innovative solution that simplifies spontaneous pneumothorax and prolonged air leaks management. P-eggy can be directly connected to any chest catheters, combining safety, comfort, and cost-effectiveness.