


AIR-LIQUID SEPARATION CHAMBER

Rapid air evacuation through a low-pressure one-way valve. Prevents liquid backflow.



ACTIVE DRAINAGE

Press to activate drainage or to verify system patency.



AIR-LEAK DETECTION

Qualitative air-leak evaluation (bubbling).

CREDITS

- [1] CTO Hospital A.O. dei Colli, Napoli, Italy
- [2] IRCCS Clinical Institute HUMANITAS Research Hospital, Milan, Italy
- [3] Cattinara Hospital, ASUGI, Trieste, Italy

Drentech WearABLE

WearABLE, drainage bag, fixing belt, adaptors
Packaging: 6 units for box

Drainage bag

Packaging: 10 units for box

Fixing belt

Packaging: 20 units for box



Via Galileo Galilei, 18 - 46025 Poggio Rusco (MN) ITALY
Tel. +39 0386 830582 | Fax +39 0386 51898
E-mail: info@redax.it | www.redax.it

Code **10829**

Code **12005**

Code **10569**

video tutorial



REV. 09-2025



Safe and wearable thoracic drainage.
Smart Design, ready-to-use.



FROM THE OR, TO HOME



Simplifies ambulatory pleural management

CRITICAL CARE



FIRST AID & AMBULANCE



No water needed, no training needed

HOME CARE



Continuous drainage is ensured in any patient position



Ideal to promote early mobilization

CARDIO-THORACIC SURGERY

